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Rev. 6 (06/03/2026)

*CERVICAL DISC DECOMPRESSION KIT WITH QMR TECHNOLOGY**

Instructions for Use
ENGLISH

CONTENT OF PACKAGING

The cervical disc decompression kit Resadisc® short, identified with product code 2615018, is composed of the following pieces:

- 1 bipolar electrode short Resadisc®;
- 1 19G introducer needle with a stylet and safety stopper;

INDICATIONS FOR USE

The Resadisc® is intended for the ablation of soft tissue in surgical procedures, such as disc decompression or treatment of protrusions.

CONTRAINDICATIONS

The use of the Resadisc® bipolar electrode is contraindicated in the following cases:

- Patients with cardiac pacemakers or with other active electronic devices;
- Pregnancy and lactation;
- Discs already performed with serious impairment of the integrity of the fibrous anulus due to its surgical removal;
- Severe disc degeneration with disc height less than 50% of the standard
- Discs subjected to metal cage implants adjacent to the one to be treated;
- Erulites of anulus with free disc fragments;
- Patients suffering from spinal infections / spinal disorders / severe osteoarthritis;
- Patients with bony deformity or severe congenital anomalies.

WARNINGS

Do not touch the distal end (tip) of the Resadisc® electrode while the Resablator generator supplies power (live electrode).

During the percutaneous surgical procedure, submit the patient to local anaesthesia, to allow for monitoring and to check if there are any signs of irritation of the nerve root.

Always check the good visibility of the Resadisc® tip via imaging before powering up.

As you perform decompression through Resadisc®, stop the procedure immediately if the patient experiences painful symptoms.

Do not extract the Resadisc® or the insertion needle from the surgical site while the Resablator generator is delivering power.

Accidental contact or immediate proximity of metal objects to the Resadisc® electrode in operation may damage the tip and cause nervous stimulation.

Electrode wear due to prolonged use may compromise the performance of the Resadisc® electrode.

PRECAUTIONS

Before using, carefully read all the documentation included in the package, including the warnings, precautions, operating instructions and operating instructions.

Telea Electronic Engineering disclaims any liability in case of improper use of the electrode.

The short Resadisc® kit is supplied in sterile packaging and is disposable. Do not resurface and / or reuse.

Cleaning, re-sterilizing or re-using the disc decompression kit may expose the patient to a serious risk of infection.

If the short Resadisc® Kit is damaged, do not use the product or try to sterilize it.

It remains at the discretion of the physician to evaluate in patients the indication of minimally invasive surgery.

Before use check that the Resadisc® electrode connector is fully inserted into the Resablator generator connector corresponding to the "BIPOLAR" output.

Avoid getting fluids in contact with the pins of the connectors.

Avoid contacting the patient's body with metallic objects.

Use or rubbing on the bony surfaces may result in excessive wear or tear of the electrode.

Always comply with fire regulations. Sparks and heat (though contained) produced in electro-surgical procedures may be potential sources of fire.

Do not use flammable agents for cleaning and disinfection of the Resablator generator.

As with other electro-surgical units, electrodes and cables can act as a path for high-frequency current. Place the cables in order to avoid contact with the patient or other conductors. Placement of other electrical equipment next to the system may cause interference.

High-frequency electro-surgical equipment may interfere negatively with the operation of other electronic equipment.

The monitoring electrodes should be positioned as far as possible from the surgical electrodes, in case of concomitant use on the same patient of high-frequency surgical equipment and physiological monitoring equipment.

It is therefore recommended to always use monitoring equipment equipped with high frequency current limiters.

We do not recommend the use of needle monitoring electrodes.

SIDE EFFECTS

Following electro-surgical intervention, though min-invasive, the tissue adjacent to the entry area could have slight iatrogenic lesions.

TRAINING REQUIREMENTS OF THE OPERATOR

These instructions and devices described in this manual must be used only by qualified medical personnel.

The operator must have acquired specific experience in general spinal surgical techniques and electro-surgery.

The user is advised to keep up-to-date on the progress made in the field of spinal electrosurgical procedures. It is also recommended to follow a training course on the short Resadisc® cervical disc decompression kit operational technique and the Resablator generator maintained by AMS Group Srl s.u.

USE MODE

Resadisc® is suitable for low-temperature decompression of the disc core in the treatment of patients with discrete hernia symptomatology and disc protrusions.

Prepare the patient for surgery according to normal preoperative procedures. Check the integrity of the packaging and its contents before use.

Prepare the generator according to the instructions provided in the Resablator's operating manual. The Resadisc® bipolar electrode is designed for operation with the only Resablator Smart generator. An adapter (code ref. 2505065) is available upon request for connection to the specific Resablator 50 generator model.

CAUTION: Before connecting the Resadisc® electrode to the generator, check the compatibility of the components of the kit: Remove the stylet and insert the Resadisc® electrode into the introducer needle by checking the smooth sliding of the stem.

The sliding between the parts may be subject to some degree of friction, therefore exerting adequate pressure without compromising the integrity of the Resadisc® electrode tip. In case of excessive friction, replace the Cervical Disc Decompression Kit with a new one.

At the end of the test, remove the Resadisc® electrode from the introducer needle and reinsert the stylet.

Connect the connector located at the proximal end of the Resadisc® bipolar electrode, at the bipolar output (labelled "BIPOLAR") of the Resablator generator.

CAUTION: The correct coupling between the electrode and the generator is done with a "click". The generator display will display screens confirming the reading of the Resadisc® electrode.

Select the bipolar mode on the Resablator front panel by pressing the button and check the activation by displaying the light at the top right of the button.

Using the "POWER" knob, adjust the output power value recommended for the method at 16.

CAUTION: Using power settings above **26** can cause undesirable effects on the tissues.

CAUTION: The resonance depth at "Quantum Molecular Resonance" of the Resadisc® electrode in the disc core is a parameter dependent on the power setting and the electrode's remaining time on the pathological tissue. Use a reference power of **16** and a dwell time of 3 seconds to achieve the desired final effect.

Pressing the pedal associated with the Resablator set with the above parameters activates the Resadisc® electrode.

In the case specific for each Resablator (Resablator 50 or Resablator SMART) model, refer to the instructions in the user manual of the generator.

OPERATING TECHNIQUE FOR CERVICAL DISC DECOMPRESSION WITH QMR TECHNOLOGY *

1. Access the disc you want by using the 19G Sterile Injector Needle to maintain site nerve structures near the target to be treated.

CAUTION: Use only the insertion needle supplied with the cervical disc decompression kit by Telea Electronic Engineering

Different needles may not have a compatible internal diameter and / or a correct length.

Telea Electronic Engineering disclaims any liability for damage resulting from the use of Resadisc® with any other needle that is not present in the package.

2. Using a fluoroscopic imaging technique, insert the needle into the patient's body and bring the needle tip over the anulus / nucleus junction, confirming the position using the anterior-posterior (AP) and lateral views. At this point, the tip of the needle lies within the disc core to decompress.

3. Using a fluoroscopic imaging technique and referring to the AP and lateral views of the cervical disc:

a. Lock the needle on the patient with the attached safety stopper. Remove the stylet and insert the Resadisc® electrode into the introducer needle.

b. Screw the Resadisc® electrode onto the introducer needle using the Luer Lock system. This ensures that the active part of the Resadisc® electrode is positioned correctly inside the disc, then

imaging verification will ensure that the tip of the Resadisc® electrode is inside the core and far enough from the vertebral plates.

c. Confirm, using fluoroscopic techniques, that the tip of the Resadisc® electrode remains within the nucleus and does not contact the annulus.

WARNING: The tip of the Resadisc® electrode must not contact the anterior annulus.

4. Once the desired position in the disc has been confirmed, perform the decompression according to the following sequence:

a. Press the Resablator footswitch, set to bipolar Quantum Molecular Resonance mode, and rotate the Resadisc® electrode 360° using the tabs above the Luer Lock system for 3 seconds.

b. To complete the disc decompression procedure, repeat step "a" two more times. Do not exceed 3 seconds of continuous pedal activation.

CAUTION: Do not dispense power while the tip of the electrode is inside the introducer needle.

CAUTION: If the patient complains of painful symptoms while activating the Resadisc® electrode, immediately stop the procedure and:

- Carefully examine the AP and side views under the fluoroscopic guide.
- Confirm that the tip of the needle is properly positioned within the rear edge of the anulus.
- Confirm the correct positioning of the Resadisc® tip inside the core.
- Refrain from proceeding until it is possible to confirm the correct positioning within the core.

If the patient complains again about the sudden onset of pain during decompression, discontinue the operating procedure.

CAUTION: If the nerve root or spinal cord comes into direct contact with the tip of the Resadisc® during decompression, serious injury to the nerve may occur.

CAUTION: If the Resadisc® movement causes contact with the vertebral body plates during decompression, severe injury to the plates may occur.

CAUTION: The bipolar electrode Resadisc® is not suitable for treating the fibrous tissue present in the anulus.

5. After the decompression of the disc core, interrupt the power supply and subsequently extract the Resadisc® electrode from the introducer needle and then the needle from the body of the patient.

WARNING: Avoid extracting the Resadisc® electrode from the trigger needle with the pedal activated.

WARNING: When the procedure is completed, in case of anchoring with fabric or excessive friction of the electrode within the 19G introducer needle, rotate the handpiece to facilitate the extraction.

WARNING: In case of excessive anchorage, if the handpiece rotation is ineffective or impossible, remove both the electrode and the 19G insertion needle.

6. Dispose of the Resadisc® kit in accordance with applicable regulations.

7. Follow the standard skin cleansing procedures and apply a sterile dressing above the needle outlet point.

8. Turn off the generator by pressing the power switch.

SAFETY CHARACTERISTICS

The safety features of the Resablator generator are described in the operating manual supplied with it.

SUPPLIES

The Resadisc® Cervical Disc Decompression Kit is sterile and is disposable. The product is only sterile if the package is sealed, untouched and without any damage. Do not use the product or try to re-sterilize or reuse it if the Resadisc® kit is damaged. Do not clean, re-sterilize or reuse as the result may be harmful or compromising on the performance of the device and expose the patient to the risk of infection.

The Resadisc® cervical disc decompression kit is sterilized by ETO.

MAINTENANCE AND INDIVIDUATION OF TROUBLESHOOTING

For instructions on how to operate and detect faults in the Resablator generator refer to the manual of the model in use.

CUSTOMER CARE

The Cervical Disc Code Decompression Kit 2615018 is guaranteed against defects in material, operation and processing for use on a single patient.

Telea Electronic Engineering it will in no way be liable, explicitly or tacitly, for any damages caused or caused by the customer or the user of the product as a result of:

- (a) Improper use, negligence and / or misuse;
- b) Repairs or modifications not made by Telea Electronic Engineering or authorized repair facilities;
- (c) Use in any mode or medical procedure other than those for which the product was designed, and any special and / or consequential special damages of any kind but related to the sale or use of the product.

This warranty replaces any other express, implied, and / or statutory warranty, including, but not limited to, any implied warranties of merchantability or fitness for a specific purpose and all other obligations or liability by Telea Electronic Engineering

COMPLAINTS

All questions regarding the quality, reliability and / or durability of this product and incident reports must be sent to Customer Service or a representative authorized by AMS Group Srl s.u. Any serious accident occurring in relation to the device must be reported to the Manufacturer and the competent authority of the Member State in which the accident occurred.

MANUFACTURER: Telea Electronic Engineering S.r.l.
Via Leonardo da Vinci 13, 36066 Sandrigo (VI)
info@teleamedical.com

LEGEND

GRAPHIC SYMBOLS USED ON LABEL			
Symbols	Symbols	Symbols	Symbols
	Do not reuse - Disposable		Expiry date
	Lot number		Sterilization of Ethylene Oxide
	Code number		Refer to the operating instructions
	Do not use if the package is not complete		Manufacturer
	Do not re-sterilize		CE certification and NB code
	Non-pyrogenic		Manufacturing date

Telea Electronic Engineering declines any liability for damages caused to the patient and / or staff if the system is used abnormally, improperly or as a result of erroneous manoeuvres and unauthorized personnel.

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HERSTELLER: Telea Electronic Engineering S.r.l.