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The screenshot shows the 'PROGRAMMI PREIMPOSTATI' screen of the i-STIM software. The interface is split into two main sections for CH1 and CH2. Each section contains a 'Programma selezionato' dropdown menu, a 'PROG:' dropdown menu, and three sliders for 'Frequenza selezionata' (F1, F2), 'Larghezza impulso selezionato' (T1, T2), and 'Tempo'. The sliders are labeled with 'Hz' or 'ms' and 'min'. A 'PROCEDI' button is located at the bottom center of the screen. The top bar includes the i-STIM logo, a home icon, and the text 'Programmi personalizzabili 8,9 e 10 Selezionare un programma su CH1 e CH2'.

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Ams Group

i-STIM
ULTIMATE BIPOLAR PENS TREATMENT

USER MANUAL

Radiofrequency generator

I-STIM

Manufactured by:

AMS Group S.r.l. s.u.

Via Europa, 12

35020 San Pietro Viminario (PD) – ITALY

VAT NO. 03638130280, Fiscal Code 03365730260

www.amsvita.com - *info@amsvita.com*


Attachments: Declaration of Compliance, Test Certificate

INDEX:

Medical Device Label:	3
Meaning of symbols:	3
Warnings:	4
Warning:	5
The use of the device is contraindicated:	5
Side effects:	5
Intended use:	6
Technical characteristics	7
ELECTRICAL SYSTEM	7
TRANSFORMER	7
MECHANICAL PROPERTIES	7
ENVIRONMENT	7
During operation:	7
During storage/transport and with the system off:	7
NEEDLES SPECIFICATIONS:	7
ENVIRONMENTAL PROTECTION - DISPOSAL	8
Package Contents:	8
Device description	9
How to use the product	10
Device installation	10
How the device works:	11
Type and selection of electrostimulation programs:	17
Maintenance and cleaning	18
Malfunctioning	18
EMC Tables (CEI EN 60601-1-2)	1
Tables 1, 2, 4, 6, 9 according to the regulations EN 60601-1-2 (2015)	1

Tables 1, 2, 4, 6, 9 according to the regulations EN 60601-1-2 (2015)..... 19

Ams Group, manufacturer of the i-stim electro-stimulator, thanks you and congratulates you for choosing our medical device. In order to ensure an effective and practical partnership, in the common interest, we would be grateful if you would like to contribute with your experience and your suggestions, to express opinions in order to constantly improve the product, where this is currently possible.

Medical Device Label:**Meaning of symbols:**

	It is compulsory to consult the instructions for use		Dispose of according to current regulations
	Manufacturer		Relative humidity
	Year of manufacture		Pressure
	Lot number		Operating temperature
	Reference code		Protection class against electrical hazards. Type BF
	Series number		CE marking and identification of the Notified Entity

Warnings:

Do not use the device until you have read and fully understood everything in this Ams Group manufacturer's user manual. Keep these operating instructions in the vicinity of the device.

Ensure the presence of the operator in the vicinity of the device during the entire time of its use. In order to avoid the risk of electric shock, this equipment must only be connected to power supply networks with a protective earth.

The device is equipped with a display/touch screen. If abnormal behavior is observed, discontinue use of the device immediately and contact technical support authorized by the manufacturer.

Keep the device out of all wet environments. Do not spill liquids of any kind on the device and do not immerse the device or any part of it in liquids. If you suspect that liquid may have entered, or is entering the device, disconnect the device from the power supply immediately. For this operation, make sure that any liquid has not entered the power cord as well, in which case, turn off the mains power switch before disconnecting it from the mains. If you suspect that liquid may have entered the device, contact technical support before using it again.

Warning:

The use of the unit described in this booklet is currently only permitted to medical personnel who are recognized and authorized to carry out medical therapeutic activities, with specific training in pain therapy, operating in public or private medical clinics or facilities and/or hospitals.

Warning: this device must only be connected to standard electrical networks.

The use of the device is contraindicated:

The PENS procedure, with i-Stim Radiofrequency Generator and i-PROBE Bipolar Needle, is contraindicated in the following cases:

- In patients with an implantable electronic device (e.g. cardiac pacemaker, defibrillator or neurostimulator).
- In patients with heart problems or epilepsy problems.
- In patients with a confirmed or presumed pregnancy.
- In patients with thrombotic problems or those undergoing anticoagulant therapy.
- In patients with infection at the site to be treated.
- In pediatric patients.

It is also recommended that the therapy be intended for patients who are alert and have mental capacity.

DO NOT apply the needle to damaged skin and DO avoid contact of the stimulation cables with damaged skin.

Side effects:

Infection, bleeding or bruising at the insertion point of the probe; Temporary increase in painful sensations.

Intended use:

The medical device known as "i-Stim" is intended to be used to interrupt the transmission of peripheral neuropathic pain by blocking the signal located in the affected nerves (analgesic effect) in association with the disposable and sterile bipolar needle "i-Probe", which can also be used to perform medicinal infusions on the target tissue.

It is an active medical device, connected to an certified medical external power supply that transforms the mains voltage from 230V ac to 5V dc. The device uses sterile disposable needles for the transmission of the electrical impulse to the patient and must be used exclusively for the treatment of human pathologies, according to the methods specified by the therapeutic protocols.

The use of the medical device described in this manual is currently only permitted for medical personnel who are recognized and authorized to perform medical therapeutic activities, with appropriate specific training in the field of pain therapy, working in outpatient clinics or medical facilities and/or public or private hospitals.

The subcutaneous electrostimulation known as Pens is a direct procedure that stimulates the peripheral nerves and nerve endings through an electrical impulse generated by the medical device and conveyed to the patient through the sterile disposable needle. The treatment is painless, repeatable, of limited invasiveness (no implanted device), can cover any anatomical area and the application on the part to be treated lasts a maximum of 30 minutes.

The stimulation is performed by means of needle probes inserted into the patient's skin. The treatment is used to relieve chronic neuropathic, nociceptive and mixed neuropathic pain, resistant to conventional therapies, particularly cases of hyperalgesia and allodynia in different areas of the body that are limited in size.

The main therapeutic indications are the following:

- Treatment of chronic peripheral neuropathy;
- Treatment of nociceptive neuropathy resistant to conventional therapies;
- Chronic postsurgical pain;
- Cervicalgia;
- Headache;
- Muscle pain (myofascial) encountered both at rest and during movement;
- Occipital / supraorbital / intercostal neuralgia
- Postherpetic neuropathy;
- Phantom limb pain.

Improper use of the device, negligence on the part of the user, failure to comply with hygiene regulations or any maneuver or use other than that expressed in this file and the protocols deposited, releases the manufacturer from all liabilities.

The set of parts that can come into contact with the patient and that, therefore, are similar to applied parts, consists of the disposable needle including the cable, the handpiece and the connector.

Position the pulse generator medical device so that it is not accessible and operable by the patient.

Use a gauze to support the cable and handpiece during treatment.

Technical characteristics

ELECTRICAL SYSTEM

- Power supply 5V dc 50/60 Hz
- Maximum current consumption 4A
- Insulation class Class 2: Type BF
- Power supply insulation Double insulation
- Transformer Medical type
- Display Tft LCD capacitive, 7-inch diagonal
- Firmware version Rev.1 - 26/06/2019

TRANSFORMER

- Brand Mean Well
- Code GSM36Bo5-P1J
- Input voltage 230V
- Output voltage 5v
- Power 22,5w
- Current 4,5A
- Type Medical

MECHANICAL PROPERTIES

- Width 285mm
- Depth 198mm
- Height 61,2mm
- Weight 1,5 kg

ENVIRONMENT

During operation:

- Temperature: 15-35°C
- Humidity: 30-75%
- Pressure: 700-1060 hPa

During storage/transport and with the system off:

- Temperature: 0-50°C for durations of less than 48 hours 0-35°C for durations longer than 48 hours
- Umidità: 5-85%
- Pressione: 700-1060 hPa

NEEDLES SPECIFICATIONS:

Bipolar needle developed with integrated infusion line; the needle is produced in various diameters of 18G, 21G, 22G and in different lengths: 35mm, 70mm, 100mm, 120mm, 200mm. The needle + generator system is able to read the impedance value (range: min. 100 Ohm, max. 1000 Ohm) of the tissue to be treated. Length of the connection cable to the medical device pulse generator is 2 metres. Needle sterilization method used: ethylene oxide.

LIST OF NEEDLE CODES:

Code/model	Description	Brand
BPP-1835IL	Bipolar Pens Probe 18G L. 35mm <i>with infusion line</i>	I-Probe
BPP-1870IL	Bipolar Pens Probe 18G L. 70mm <i>with infusion line</i>	
BPP-2135IL	Bipolar Pens Probe 21G L. 35mm <i>with infusion line</i>	
BPP-2170IL	Bipolar Pens Probe 21G L. 70mm <i>with infusion line</i>	
BPP-21100IL	Bipolar Pens Probe 21G L. 100mm <i>with infusion line</i>	
BPP-21120IL	Bipolar Pens Probe 21G L. 120mm <i>with infusion line</i>	
BPP-2235IL	Bipolar Pens Probe 22G L. 35mm <i>with infusion line</i>	
BPP-2270IL	Bipolar Pens Probe 22G L. 70mm <i>with infusion line</i>	
BPP-22100IL	Bipolar Pens Probe 22G L. 100mm <i>with infusion line</i>	
BPP-22120IL	Bipolar Pens Probe 22G L. 120mm <i>with infusion line</i>	
BPP-2235	Bipolar Pens Probe 22G L. 35mm	
BPP-2270	Bipolar Pens Probe 21G L. 70mm	
BPP22100	Bipolar Pens Probe 22G L. 100 mm	
BPP-22120	Bipolar Pens Probe 22G L. 120 mm	
BPP-22200	Bipolar Pens Probe 22G L. 200mm	

ENVIRONMENTAL PROTECTION - DISPOSAL

The packaging of the device must be disposed of by waste recycling. Metal parts must be disposed of through the metal disposal service. Plastics, electrical components, and electronic boards must be disposed of as electronic waste. Disposal must be carried out in accordance with national directives and laws using appropriate waste disposal companies. All supplied accessories must also be disposed of in accordance with the regulations.


Warning: All parts of the device must be disposed of to avoid the risk of environmental pollution, do not disperse in the environment.

Package Contents:

- N°1 I-stim Medical Device
- N°1 Sterile single-use bipolar needle with connection cable, model: MDL (see needle code list page 8)
- N°1 Medical power supply model: GSM36B05-P1J Mean Well
- N°1 User Manual


Device description

Figure 1: Front panel

- 1. Touchscreen display
- 2. Needle input N°1
- 3. Manual pause button
- 4. Needle input N°2

Figure 2: Back panel

- 5. Usb ports (service only)
- 6. Power supply connection
- 7. Main switch ON/OFF

How to use the product

This device is designed to generate electrical impulses that are delivered to the patient via sterile single-use needles. It needs to be powered with 5V electric current through an external medical power supply (model: GSM36B05-P1J Mean Well) that transforms the mains voltage from 230V to 5V dc, plugged into the appropriate socket. **The electrical system must be in compliance with regulations.** The power supply is part of the medical device, therefore the use of power supplies other than those provided by the manufacturer may cause an increase in emissions or a decrease in the integrity of the electromedical device. It may also damage the device and void the warranty. Therefore, use the medical device only with the medical power supply model: GSM36B05-P1J Mean Well.

The medical device only works with MDL branded needles supplied by the manufacturer. Electrodes or needles other than those supplied by the manufacturer may damage the device and void the warranty.

Device installation

The installation of the device is a simple operation that can be performed by both the end user and the distributor.

Remove the medical device from its box by removing the packaging. Check the contents of the package (see page 8 "package contents").

Steps to be taken to install the device:

1. After removing the device from its packaging, place it on a smooth, horizontal surface near a power outlet. The back of the device must be free of objects that may obstruct the connections.
2. Take the power supply supplied with the device and connect it on one side to the device and on the other side to the 230V power outlet.

IMPORTANT:

Do not place the device on other equipment or vice versa. Position the device so that the connecting devices (needles, power supply) cannot disconnect from it, interrupting treatment.

In case of emergency, press the manual stop button or switch off the power by pulling out the power supply connector on the rear side of the device. These operations interrupt the operation of the device.

Do not apply needles near the heart, eyes, head and any sensitive areas that may cause cardiac fibrillation in the patient.

ALWAYS RETURN THE DEVICE TO ST/BY AFTER EACH TREATMENT

THE MANUFACTURER DECLINES ALL RESPONSIBILITY FOR THE IMPROPER USE OF THE DEVICE AND FOR ANY DAMAGE OR INJURY CAUSED BY NEGLIGENCE AND/OR INEXPERIENCE OF THE OPERATOR.

DO NOT PLACE THE MEDICAL DEVICE IN SUCH A WAY THAT IT IS DIFFICULT TO DISCONNECT THE DEVICE FROM THE POWER SUPPLY.

How the device works:


1. Insert the plug of the supplied power supply unit into a mains socket correctly connected to the earth; insert the connector of the cable itself into the appropriate socket on the machine.
2. Activate the device by pressing the ON/OFF button on the back of the device.
3. Verify that the display turns on and that it runs a control sequence for approximately 20 seconds. Screen 1 then appears.
4. You can stop stimulation at any time by pressing the manual stop button on the front panel.

Screen 1: When the device is switched on after the system control sequence, the Amsgroup logo appears for a few seconds.

Screen 1

Screenshot 2: The i-stim logo is displayed and the device performs a self-check illustrated by the bar loading and the animation in the lower left corner.

Screen 2

Screen 3: The home page of the device is displayed where you can access the program list and channel screen. Only from this window you can turn off the device **IN A SAFE WAY** by using the button with the power off symbol on the top bar.

Screen 3

Screen 4: it appears when you press the "PROGRAMS LIST" button from the home page. The programs list is displayed, pressing the buttons related to the program number, the detailed description of the program appears on the right of the window. There are 7 preset programs and 3 customizable ones. By pressing the "go to preset programs or go to customizable programs" keys, you can access the channels screen. From this page you can return to the home page by pressing the home icon on the top bar.

Screen 4

Screen 5: the screen of the preset programs applicable to CH1 and CH2 (channel1 and channel2) is displayed.

From the drop-down menu you can also choose two different programs, one for each channel, and use them simultaneously. Each program has specific characteristics described in the table on page 15.

To continue to the next screen to start treatment, press the Proceed button, but first connect the needle and select a program. Without these conditions, it is not possible to continue further. A warning will appear to select an electrostimulation program and the connected needle check will be displayed on the screen.

By pressing the customizable programs button on the top bar you can switch to the customizable programs screen without returning to the program list page.

From this page you can return to the home page by pressing the home icon on the top bar.

From this page you can return one page back by pressing the arrow icon on the top bar.

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PROCEDI

Screen 5

Screenshot 6: Displays the screen of customizable programs applicable to CH1 and CH2 (channel1 and channel2). You can adjust the frequency, pulse amplitude and time on both phase 1 and phase 2.

From the drop-down menu you can also choose two different programs, one for each channel, and use them simultaneously. Each program has specific characteristics described in the table on page 17.

In order to continue to the next screen to start the treatment, press the Proceed button, but first you have to connect the needle and select a program, without these conditions you cannot continue any further. A warning will appear to select an electrostimulation program and the check of the connected needle will be displayed on the screen.

Pressing the Preset Programs button on the top bar allows you to switch to the Preset Programs screen without returning to the program list page.

From this page you can return to the home page by pressing the home icon on the top bar.

The screenshot displays the 'PROGRAMMI PREIMPOSTATI' screen of the i-STIM software. The interface is split into two main sections for CH1 and CH2. Each section includes a 'Programma selezionato' dropdown menu, a 'SELEZIONE PROGRAMMA' section with a 'PROG:' dropdown, and three sliders for 'Frequenza selezionata' (F1, F2), 'Larghezza impulso selezionato' (T1, T2), and 'Tempo'. The units for frequency are Hz, for pulse width are ms, and for time are min. A 'PROCEDI' button is located at the bottom center of the screen.

Screen 6

Screen 7: displays the start screen of the previously selected electrostimulation program on both CH1 and CH2 (channel 1 and channel 2), even at the same time. Before activating the stimulation, the user must validate the needle by pressing the button on the top bar "press here to validate the needle".

The user can adjust the operating voltage from 0V to 5V, according to the therapy to be performed and can visualize in real time the impedance value generated by the positioning of the needle in the patient's skin thereby starting therapy when the impedance value meets the parameters defined by the user.

Press the start button to start stimulation. If the needle is not connected, stimulation does not start and a warning appears. When the stimulation starts, the time is marked by an animated bar and the treatment minutes are counted.

During stimulation you can change the voltage from 0 to 5 Volts in real time. During the stimulation it is possible to visualize in real time: the set program, the current, the impedance.

To stop stimulation press the Stop button.

From this page you can return to the home page by pressing the home icon on the top bar.

From this page you can return one page back by pressing the arrow icon on the top bar.

Screen 7

WARNINGS:

Warning 1: Pressing the "Proceed" button, from the program selection screen without selecting any program, the following warning appears: select a program for channel 1 or channel 2. You will not be able to proceed further if you do not select a program for each channel.

Warning 2: If you press the "start" button from the "stimulation activation" screen, without connecting the needle, the following pop up appears: connect and/or validate a needle on channel 1 or channel 2.

Warning 3: Press here to validate the needle, if the needle is new the following pop up appears: needle on channel 1 or channel 2 connected and ready to use. If the needle has already been used, the following pop up appears: needle not connected correctly or already used on channel 1 or channel 2.

The following pop up appears during stimulation: stimulation in progress.

Note 4: If you press the "home" button or the "back" button during stimulation, the following pop up appears: stop stimulation.

Warning 5: The following warning appears at the end of stimulation: end of stimulation on Ch1 or Ch2.

Warning 6: When you press the manual pause button the following pop up appears: Turn off the device? Pressing "yes" will turn off the device, pressing "no" will return to the home page.

Warning 7: pressing the power off button the following pop up appears: Stimulation on channels 1 and 2 paused; to resume treatment press ok.

WARNINGS:

The medical device covered by this user manual must be used by medical personnel who are experienced in pain therapy.

Use only the power supply and needle provided by the manufacturer. Therapy should only be carried out on adult persons and attention should be paid to contraindications.

In case of side effects in the patient, discontinue therapy temporarily.

Do not reuse needles or use them for other purposes. Only connect the device to standard electrical installations.

Read the user manual carefully before using the device. Dispose of in accordance with applicable regulations.

If the needle is not connected, stimulation will not start. If you do not stop stimulation by pressing STOP, you cannot return to the homepage or go back one screen.

Type and selection of electrostimulation programs:

The i-stim medical device is designed to work with two needles connected at the same time because it has two separate and independent channels. It is therefore possible to act on two different areas of pain in the patient.

The manufacturer has prepared 7 predefined electrostimulation programs and 3 programs where the parameters of applied voltage and pulse duration can be customized by the user.

The maximum output values from needles are: 5V voltage, 120Hz frequency, 2ms pulse amplitude, 60sec single frequency time, 30min treatment duration, 50mA current.

Below is a table showing the parameters of the electrostimulation programs defined by the manufacturer:

I-STIM PROGRAM LIST					
Prog. N°	Frequency	Pulse width	Single Frequency Time	Treatment timing	Note
P1 (std)	2Hz	1ms	NA	30min	Fixed Frequency, Pulse width, and treatment time
P2 (Std)	100Hz	0,2ms	NA	30min	
P3	2Hz	1ms	3s	30min	Alternate stimulation: alternation of 2phases with different frequency and pulse width
	100Hz	0,21ms	3s		
P4	F1 - 2Hz	F1 - 1ms	2Hz pulse reach 100Hz pulse in 60sec	30min	Frequency Modulation: Frequency is automatically modified between min-max range and pulse width is adapted to the frequency
	F2 - 100Hz	F2 - 0,21ms			
P5	2Hz	Min. 0,21ms	60 sec	30min	Pulse width modulation: fixed frequency but pulse width change
		Max. 1ms			
P6	100Hz	Min. 0,21ms	60sec	30min	Low Frequency
		Max. 1msec			
P7	100Hz	0,21ms	0,5sec with and consequently pause of 0,5ms	30min	High Frequency Peak
P8	F1 - 2Hz	F1 - 1ms	See note	30min	Alternate Stimulation: 2 phases with different frequencies and pulse width follow each other in few seconds. Customized parameters: Phase1: from 1 to 120Hz with pulse width from 0,2ms to 2ms; Phase2: from 1 to 120Hz with pulse width from 0,2ms to 2ms;
	F2 - 100Hz	F2 - 0,21ms			
P9	Set from 1 to 120Hz	Set from 0,2 and 2ms	0,5sec e con pausa a seguire da 0,5sec	30min	Custom
P10	Set from 1 to 120Hz	Set from 0,2 and 2ms	0,5sec e con pausa a seguire da 0,5sec	30min	Custom

Maintenance and cleaning

Maintenance operations of the medical device for the operator are limited to external cleaning.

Pay particular attention to the touchscreen display; to avoid scratches you can clean it with a commercially available microfiber cloth; do not use solvents or abrasive cleaners.

Prevent liquids from spilling on the device to prevent them from getting inside. If this happens accidentally, turn off the device and contact the technical support center immediately.

Disposable needles are supplied in sterile packaging. Discard the needle after use in special containers for disposal, in accordance with the applicable regulations. Needles are not reusable or re-sterilizable.

The manufacturer does not provide for routine maintenance of the I-stim medical device, in case of failure contact the service center described below.

Authorized service centers:

ECO3 s.n.c.

Via Ceresole, 1 10032 Brandizzo TO, Italy

Tel. +39 0119139673 e-mail: eco3ozono@gmail.com

Warning!

Improper use of the device or tampering with it by personnel who are not authorized by the manufacturer will invalidate the warranty.

When shipping the device, you should use the original packaging and contact the service center for information on how to send it.

Malfunctioning

What to do before contacting the technical support center

IF IT DOES NOT TURN ON OR DOES NOT RESPOND TO COMMANDS

- Turn the device off using the rear main switch, wait at least 30 seconds and then turn it on again.
- Repeat the operation if necessary.

If problems persist, please contact support.

EMC Tables (CEI EN 60601-1-2)

Tables 1, 2, 4, 6, 9 according to the regulations EN 60601-1-2 (2015)

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions		
The i-stim device is intended to operate in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment
CISPR 11 RF emissions	Group 1	The i-stim medical device uses RF energy only for its internal operation. As a result, its RF emissions are very low and are unlikely to cause any interference in nearby electronic devices.
CISPR 11 RF emissions	Class B	The medical device is suitable for use in environments such as medical clinics and those directly connected to a public low-voltage power supply that supplies buildings used for domestic purposes.
Harmonic emission	Class A	
Voltage/Flicker Emissions fluctuations of	Compliant	

Table 2

Immunity Test	Test level IEC 60601	Level of compliance	Electromagnetic environment
Electrostatic discharge	± 8 kV on contact ± max 15 kV in air	± 8 kV on contact ± max 15 kV in air	Floors must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Transient / electrical pulse sequence	± 2 kV for power supply lines 100 kHz repetition frequency	± 2 kV common mode 100 kHz repeat rate	The quality of the mains voltage should be that of a typical commercial or hospital environment.
Surges	± 2 kV phase (i)-ground ± 1 kV phase-(s)	± 2 kV common mode ± 1 kV differential mode	The quality of the mains voltage should be that of a typical commercial or hospital environment.

Voltage drops, short interruptions and voltage variations on power supply input lines	< 0.5% Ut (95% drop) for 0.5 cycles 40% Ut (60% drop) for 5 cycles 70% Ut (30% drop) for 25 cycles < 5% Ut (95% drop) for 5 seconds	< 5% Ut (95% drop) for 0.5 cycles 40% Ut (60% drop) for 5 cycles 70% Ut (30% drop) for 25 cycles < 5% Ut (95% drop) for 5 seconds	The quality of the mains voltage should be that of a typical commercial or hospital environment. There are no problems in the event of a mains voltage interruption.
Magnetic field at mains frequency (50/60 Hz)	30 A/m	30 A/m	Network frequency magnetic fields must have levels that are typical of a location in a commercial or hospital environment.

Table 4

Immunity test	Test level IEC 60601	Level of compliance	Electromagnetic environment (distance of a transmitter from the equipment)
Conducted RF	3 Veff	3 Veff	$d = 1,17 \cdot \sqrt{P}$ where P is the maximum rated output power of the transmitter in W according to the transmitter manufacturer and d is the separation distance in metres
Irradiated RF	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	from 80 MHz to 800 MHz from 800 MHz to 1 GHz from 1 GHz to 2,7 GHz where P is the maximum rated output power of the transmitter in W according to the transmitter manufacturer and is the separation distance in meters

Table 6

The device is intended to operate in an environment where radiated radio frequency interference is controlled. Health care professionals using the device can help prevent electromagnetic interference by ensuring a minimum distance between mobile and portable RF communication equipment (transmitters) and the device as recommended below, in relation to the maximum output power of radio communication equipment.

Maximum output power of the specified Transmitter (W)	Separation distance to transmitter frequency (m)			
	From 15 kHz to 80MHz	From 80MHz to 800MHz	From 800MHz to 1GHz	From 1GHz to 2,7GHz
0,01	0,12	0,04	0,07	0,23
0,1	0,37	0,11	0,22	0,74
1	1,17	0,35	0,70	2,33
10	3,69	1,11	2,21	7,38
100	11,67	3,50	7,00	23,33

Table 9

Immunity to wireless RF communication devices.

Test frequency (MHz)	Band width (MHz)	Service	Modulation	Maximum applied power (W)	Distance (m)	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulsed modulation on 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz 1 kHz	2	0,3	28
710 745 780	704 - 787	LTE Band 13, 17	Pulsed modulation on 217 Hz	0,2	0,3	9
810 870 930	800 - 960	GSM 800/900 TETRA 800, IDEN 820, CDMA 850, LTE BAND 5	Pulsed modulation on 18 Hz	2	0,3	28
1720 1845 1970	1700 - 1990	GSM 1800, CDMA 1900, GSM 1900 DECT, LTE BAND 1,3,4,25, UMTS	Pulsed modulation on 217 Hz	2	0,3	28
2450	2400 - 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE BAND 7	Pulsed modulation on 217 Hz	2	0,3	28
5240 5500 5785	5100 - 5800	WLAN 802.11 a/n	Pulsed modulation on 217 Hz	0,2	0,3	9

The minimum separation distances for the highest immunity test levels are calculated using the following equation:

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, E is the immunity test level expressed in V/m.

Warranty conditions

The medical device is covered by the warranty for 12 months from the date of purchase certified by a tax receipt, invoice or transport document showing the reference data of the device and the serial number. The manufacturer does not recognize other forms of warranty issued by dealers, resellers or third parties. In the absence of the above-mentioned documents proving the date of purchase, any intervention is at the total expense of the purchaser.

- Warranty means the repair or replacement free of charge of parts or components recognized as defective due to manufacturing or material defects of the unit recognized by Ams Group.

- The device will be repaired in one of the laboratories authorized by the manufacturer or at the manufacturer's premises. The costs and risks of transport to and from the abovementioned laboratories shall be borne in full by the purchaser.

- If a "call out" repair is required, this will only be carried out if technically possible, the applicable fixed call charge and any transport costs must be paid if the repair requires the device to be returned to the laboratory.

- This warranty operates only if the device has been purchased and used in one of the member countries of the European Union with the exception of the so-called "free zones" and if it has been used in accordance with the instruction manual and the warning labels located on it.

- The warranty does not apply to damage caused by carelessness, misuse, mistreatment, deterioration, accidents, vandalism, natural phenomena, nor to damage due to installation, adaptation or modification that has not been authorized by personnel recognized by the manufacturer and/or authorized from time to time by the same manufacturer. Furthermore, it does not apply to damage caused by incorrect use or contrary to the technical and/or safety measures required in the country in which the device covered by this warranty is used or in any case, as provided for in the instruction manual.

- As soon as the purchaser discovers a defect in the device, he must promptly contact the local dealer or branch (where applicable) or directly the manufacturer or the national distributor in the EU country where the warranty is invoked.

- The sale of the device by third parties without the authorization of the manufacturer in countries where there are different regulations concerning electrical safety certification or which require documentation other than that attached in accordance with EU legislation, releases the manufacturer from any liability.

i-STIM
ULTIMATE BIPOLAR PENS TREATMENT

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The screenshot shows the 'PROGRAMMI PERSONALIZZABILI' screen in the i-STIM software. The interface is split into two main sections for CH1 and CH2. Each section contains a 'Programma selezionato' dropdown menu and a 'prog' button. Below these are four rows of parameters: 'Frequenza' (F1 Hz, F2 Hz), 'Larghezza impulso' (T1 ms, T2 ms), 'Tempo singola frequenza' (T), and 'Durata trattamento' (T min). A 'Descrizione' label is at the bottom of each column. At the bottom center is a large 'PROCEDI' button.

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