BISTURI A RISONANZA QUANTICA MOLECOLARE
QUANTUM MOLECULAR RESONANCE GENERATOR

RESABLATOR 50

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ENG *These operating instructions are an integral and fundamental part of the VESALIUS® device, because they describe its operation and use, therefore they must be carefully read prior to proceeding with the installation and use of the device itself.*
All the safety instructions or notes must be observed.
This manual and the equipment described therein must only be used by qualified personnel and physicians, specifically trained for the particular surgical techniques and procedures to be applied.
Further technical information can be found in the VESALIUS®Service Manual.
No part of this document may be copied, reproduced or translated in another language without the written consent of TELEA.

Caution: Federal law restricts this device to sale by or on the order of a physician. 21CFR801.109 (b)

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www.teleamedical.com

INDEX

1	SAFETY.....	35
1.1	DEFINITIONS.....	35
1.2	WARNINGS AND OPERATING INSTRUCTIONS	35
1.2.1	GENERAL SAFETY INFORMATION.....	35
1.2.2	EXPLOSION OR FIRE HAZARDS	35
1.2.3	ACTIVE ACCESSORIES	35
1.2.4	HAZARDOUS SITUATIONS IN RESABLATION MODE	36
1.2.5	HAZARDOUS SITUATIONS IN NERVE STIMULATION MODE	36
1.2.6	HAZARDS ASSOCIATED WITH THE NEUTRAL PLATES	36
1.3	MEDICAL AND SURGICAL INFORMATION.....	37
2	PRESENTATION.....	38
2.1	MODEL TREATED	38
2.2	THE RESABLATOR50	38
2.2.1	TECHNOLOGICAL PRINCIPLES OF THE VESALIUS® RESABLATOR50	39
2.2.2	CHARACTERISTICS OF THE RESABLATOR50.....	39
2.3	TECHNIQUES AND FUNCTIONS USED	39
2.4	THE PROFILE OF THE OPERATOR AND THE ENVIRONMENT WHERE IT IS USED	40
3	DESCRIPTION OF THE DEVICE.....	41
3.1	SYMBOLS	41
3.2	FRONT PANELS DESCRIPTION	42
3.3	DESCRIPTION OF THE FRONT PANEL CONTROLS	42
3.3.1	SELECTOR SWITCHES AND BUTTONS SECTION	42
3.3.2	OUTLETS SECTION	43
3.3.3	FUNCTIONS SETTING SECTION	44
3.3.4	THE CONTROL SECTION.....	45
3.3.5	THE DISPLAY SECTION	45
3.4	REAR PANELS DESCRIPTION	46
3.4.1	SERIAL NUMBER.....	46
3.4.2	CUSTOMIZATION.....	47
3.5	DESCRIPTION OF THE REAR PANEL CONTROLS.....	47
3.6	STANDARD EQUIPMENT OF THE DEVICES	48
4	TECHNICAL SPECIFICATIONS.....	49
4.1	INPUT TECHNICAL DATA.....	49
4.1.1	FUSES	49
4.2	OUTPUT TECHNICAL DATA	49
4.2.1	OUTPUT POWER	49
4.2.2	OUTPUT VOLTAGE	49
4.3	OUTPUT CHARTS	50
4.3.1	FEATURES OF THE OUTPUT POWER.....	50
4.3.2	FEATURES OF THE OUTPUT CURRENT AND VOLTAGE	50
4.4	DEVICE FEATURES	51
4.4.1	SELF-TEST AT START-UP	51
4.4.2	ALARMS.....	52
4.4.3	PROTECTION FROM DEFIBRILLATION	52

4.5	NERVE STIMULATOR TECHNICAL DESCRIPTION	52
4.5.1	GENERAL INFORMATION	52
4.5.2	WAVEFORM	53
4.5.3	SAFETY AND PRECISION	53
4.6	CLASSIFICATION	54
4.7	INTENDED USE	54
4.8	ELECTROMAGNETIC COMPATIBILITY	54
4.9	FEATURES OF ACCESSORIES	56
5	INSTALLATION, USE AND WARNINGS	57
5.1	CONDITIONS OF USE, TRANSPORT AND STORAGE	57
5.2	INSTALLATION OF THE DEVICE	57
5.3	CONNECTION OF ACCESSORIES.....	57
5.4	GENERAL INSTRUCTIONS OF USE FOR THE NEUTRAL PLATES.....	57
5.4.1	DISPOSABLE PLATE.....	57
5.4.2	ELASTIC BAND	58
5.5	PREPARATION FOR USE	58
5.5.1	USE AND START-UP PROCEDURE	58
5.5.2	SHUTDOWN PROCEDURE.....	58
5.6	SAFETY PRECAUTIONS.....	59
6	CLEANING, DISINFECTION AND STERILISATION	60
6.1	DEVICE CLEANING	60
6.2	CLEANING AND STERILISATION OF THE ACCESSORIES.....	60
6.2.1	CLEANING.....	60
6.2.2	STERILISATION.....	60
6.2.3	STORAGE	60
7	MAINTENANCE	61
8	TROUBLESHOOTING	61
9	WARRANTY	62

1 SAFETY

1.1 DEFINITIONS

For the purposes of this manual, the following terminology will be used:

WARNING: indicates a potentially hazardous situation, which, unless avoided, might result in death or serious injury of the patient or the operator.

CAUTION: indicates a potentially hazardous situation, which, unless avoided, might result in slight injury of the patient or the operator.

PRE-CAUTION: indicates a potentially hazardous situation, which, unless avoided, might result in damage to property.

NOTE: refers to an indication given by the manufacturer, considered important or useful for the use of the RESABLATOR50 device. This information is not correlated to a risk or a dangerous situation.

1.2 WARNINGS AND OPERATING INSTRUCTIONS

The use of RESABLATOR50 must only be assigned to qualified medical staff, specialized and adequately trained to use it, as described in this manual. Said personnel must also have the appropriate and necessary knowledge of anatomy and neurophysiology and must have read and fully understood this manual, as well as any protocol in use for performing peripheral nerve blocks. Read the instructions carefully before use.

1.2.1 GENERAL SAFETY INFORMATION

6. The power output and the stimulus voltage/current must be as low as possible for the intended objectives;
7. The patient must not come into contact with metal parts that are connected to ground or have a high coupling capacitance to ground (e.g. operating tables, metal supports, etc.).
8. **CAUTION:** To avoid the hazard of electric shock, RESABLATOR50 must only be connected to power supply networks with protective earth system.
9. If the RESABLATOR50 is not used for a certain period of time, or if you do not know the procedure for selecting the voltage to be delivered, it is advisable to start with a low voltage, increasing gradually until you reach the intensity that allows you to achieve the desired result at optimum level, obviously following the literature on this subject.
10. A malfunction of the RESABLATOR50 device could result in an unwanted increase of the output power and of the stimulus voltage/current.

1.2.2 EXPLOSION OR FIRE HAZARDS

39. The use of flammable anaesthetics, or oxidizing gases such as nitrous oxide (N₂O) should be avoided in case of thoracic or brain interventions, unless they can be sucked with a special tool. Flammable substances used for cleaning purposes, or as disinfectants, or as solvents of adhesives should be allowed to evaporate before intervening.
40. There is a risk of flammable solutions stagnating under the patient or in cavities of the human body. The fluids deposited in these areas must be dried or removed prior to using the device. There is an obvious risk of ignition of the endogenous gas. Some materials, like cotton wool or gauze, when impregnated with oxygen, can catch fire due to sparks produced by the RESABLATOR50 in normal operating condition.
41. In endoscopic procedures, check the patient's internal gas before using any of the electrodes of RESABLATOR50.
42. The instruments of the RESABLATOR50 should not be used in the presence of solid, liquid or gaseous combustible or explosive substances.

1.2.3 ACTIVE ACCESSORIES

43. Replace the electrodes if the isolating protective cover is damaged. If the metal is exposed there is a risk of current leakage that can injure the patient.
44. Replace the accessories when they or their cables show signs of deterioration.
45. Prior to proceeding with the connection and replacement of the electrodes and handles, check that the machine is turned off.
46. The connecting cables with the electrodes/accessories of the RESABLATOR50 and monitoring appliances should be positioned, as far as possible, so they do not touch the patient or other cables.

47. The RESABLATOR50 uses both monopolar and bipolar modes of operation. The outputs are available on separate connectors, therefore they cannot be mixed up thanks to the different shapes and to the number of contacts of the connectors. Pay particular attention when connecting the parts installed, to avoid damaging the panel connectors. If the connectors are mixed up on purpose, the front panel will be damaged.
48. Avoid placing inactive accessories on or near the patient, but store them in a secluded place and away from the patient.
49. Do not keep two handles connected to the device at the same time, even if they are inactive.
50. In case of surgery on body parts having a relatively small section, it is recommended to use bipolar techniques, so as to avoid damaging other portions of tissue.
51. **The operating features are guaranteed only when using original accessories.**

1.2.4 HAZARDOUS SITUATIONS IN RESABLATION MODE

52. The RESABLATOR50 and the monitoring equipments should not be connected to the same fixed output power supply (use independent cables for the connection to the power supply).
53. Avoid skin-to-skin contact by using dry surgical gauze.
54. Use monitoring systems which incorporate limiting devices for high frequency current.
55. **Do not operate under any circumstances a second plate or neutral electrode (including models of Megadyne plates), even if it is connected to other equipment (especially if they are of type BF) that is switched off.**
56. When used on the same patient simultaneously, the RESABLATOR50 and physiological monitoring equipment, all monitoring electrodes that do not contain protection resistors or high frequency inductors should be placed as far as possible from the electrodes. It is also recommended to position the + and - monitoring electrodes as close as possible. Needle monitoring electrodes are not recommended.
57. There is a possible risk for patients carrying a pacemaker or pacing electrodes, because there can be an interference with the action of the pacemaker, or the pacemaker itself can be damaged. In case you have doubts read the information given in paragraph 1.3 and ask the Cardiology Section for advice.
58. The device may interfere with other electronic and/or medical equipment that is old or non compliant with the regulations in force at the time when the RESABLATOR50 is operating.

1.2.5 HAZARDOUS SITUATIONS IN NERVE STIMULATION MODE

59. The use of the device in **NERVE STIMULATION** mode may interfere with the operation of pacemakers or other stimulators implanted or anyway connected to the patient. So do not use the stimulator without first consulting a specialist in these devices.
60. Do not use **NERVE STIMULATION** mode in the presence of flammable gases.
61. Check that the cables do not come into contact with other cables or directly with the patient.
62. Avoid operating near shortwave or microwave therapy devices.
63. Do not use **NERVE STIMULATION** mode in case of malfunction of the RESABLATOR50, or in case liquids are spilled on the device. Dry carefully after cleaning.
64. Do not use non original accessories. The use of non original accessories may change the emission features of the electrical signal, consequently altering the safety features of the function itself.
65. **The use of NERVE STIMULATION mode in the chest area is not recommended, due to the high possibility of creating cardiac fibrillation, because the electrodes are applied on the heart or close to it.**
66. Simultaneous connection of the patient to the disposable neutral plate (**RESABLATION** mode) and the conductive bracelet (**NERVE STIMULATION** mode) may cause scalding at the said bracelet.
67. We insist again on using always original accessories.

1.2.6 HAZARDS ASSOCIATED WITH THE NEUTRAL PLATES

68. If the neutral plate is used:
 - j. see the information given in paragraph 5.4;

- k. do not use the insulated neutral plate unless you can position it correctly. In this case, avoid the monopolar technique and use a bipolar technique instead;
- l. the neutral plate must be connected in a reliable manner with the patient's body, possibly on its entire surface and as close as possible to the intervention point;
- m. make sure the neutral plate is correctly positioned. As far as possible, avoid placing it on scars, fat areas, skin folds, hairy areas, bony prominences, metal implants or close to other cables or electrodes. Pay particular attention in case you treat children or thin patients;
- n. Do not allow the cable to wind and come in contact with metal parts or the patient. Avoid coiling the cables and do not, under any circumstances, wrap them around the plate or the patient;
- o. No disinfectant, rinsing liquids, physiologic solutions or other current dispersing elements should penetrate between plate and the patient's body. Also, you need to make sure that the neutral plate is correctly positioned;
- p. to avoid any interference, keep the heart monitoring electrodes as far as possible from the neutral plate and from the intervention area.
- q. a noticeable decrease in the output power or a malfunction of the RESABLATOR50 when it is set for normal operation may indicate a faulty application of the neutral plate, or a bad contact in the connections of the same;
- r. If the patient is repositioned during the procedure, check again the contact of the plate and all the connections.

69. Do not cut, round off or modify in any way the neutral plate.

70. Although the RESABLATOR50 is equipped with measures aimed at reducing the neuromuscular stimulation, such possibility cannot be completely eliminated when electric arcs are produced. The user is aware that, in sensitive structures, neuromuscular stimulation can still happen, involving a secondary risk caused by muscular contractions.

71. The RESABLATOR50 use a monitor for electrode continuity and not a monitor for the quality of the contact. Therefore, any loss of contact between the neutral electrode and the patient will not trigger any alarm tone. Instead, a plate detachment alarm will be triggered in case the connection with the neutral electrode is interrupted, when the device is operating in **RESABLATION** mode .

Do not use the RESABLATOR50 simultaneously with other electrical HF surgical units on the same patient.

1.3 MEDICAL AND SURGICAL INFORMATION

Use the RESABLATOR50 with extreme caution in patients with internal or external pacemakers. The use of the RESABLATOR50 unit may interfere with the operation of pacemakers or other stimulators implanted or anyway connected to the patient.

Therefore, do not use the RESABLATOR50 device without the favorable opinion of a specialist physician, because the interference generated during operation may cause the pacemaker to function in an asynchronous manner, or may block its operation entirely. Ask the manufacturer of the pacemaker, or the cardiology section of the hospital for further information in case you have to use the RESABLATOR50 on patients carrying cardiac pacemakers.

It is recommended to:

- never intervene in the chest area (especially close to the pacemaker or the heart).
- avoid using radiosurgery on unstable patients, highly dependant on the pacemaker;
- check the peripheral pulses before, during and after the intervention;
- maintain the duration of each emission under 5-10 seconds;

If during the intervention the power output is below the usual, do not increase the high frequency output power without having previously performed the following checks:

- check that the neutral plate is correctly positioned;
- check the handle and the connection with the electrode;
- check that the foot switch is operating properly.
- check that the cable insulation is intact.

2 PRESENTATION

2.1 MODEL TREATED

This manual refers to RESABLATOR50, which is a model of the Vesalius® line:

CODE	DESCRIPTION
2501011	RESABLATOR50

This new model RESABLATOR50 incorporates a dual function, that is, the usual **RESABLATION** mode, which allows the device to be used as a Quantum Molecular Resonance generator, is backed-up by **NERVE STIMULATION** mode, which allows the RESABLATOR50 to act as a neuro-sensory stimulator to the purpose of detecting the presence of nerves prior to proceeding with the electrosurgical procedure.

So the extreme versatility of the new RESABLATOR50 combines the features of two different devices, being able to proceed first with the stimulation and then with the surgical procedure using the same parts applied (active accessories) for both modes of operation.

2.2 THE RESABLATOR50

The name *QUANTUM MOLECULAR RESONANCE* is given by the particular manner in which the energy delivered by the device interacts with the biological tissue.

Any form of energy is transmitted from a source to a user in the form of "packets" called *energy quanta* and the value of the energy transported by the quanta depends on the frequency of the source that has produced them. In its turn, the reaction that the incident energy causes on the targeted tissue has to do with the value of the energy associated with each quantum. When the energy quanta hit the recipient, the same undergoes the following type of action:

- if the quantum energy value is different from the binding energy of molecule hit, you get only an increase of the kinetic energy of the molecule, without breaking their bond, which results in an increase in temperature;
- if the energy value of the quantum is equal to the binding energy of the molecule affected, then the entire energy of the quantum is used to break the above mentioned bond, without increasing the kinetic energy and therefore without increasing the temperature.

With the RESABLATOR50, by selecting the cut function, the energy quanta transmitted reach the tissue having the same value as the binding energy of the molecules the tissue consists of. That way, the incident quanta energy causes the bonds of the tissue molecules to resonate, neatly breaking them in a sharp way. Consequently, the rise of the tissue temperature is minimum, because the energy is mainly transferred as potential rather than kinetic.

Molecular bond cleavage leads, in its turn, to the rupture of the cell membranes (cutting effect), significantly reducing the thermal damage. The cut obtained is precise and very delicate, leaving the surrounding tissue perfectly intact.

Unlike the more traditional electric or radio surgery scalpels, whose operating principle is a mere transfer of thermal energy (heat generated by the passage of electric current), the RESABLATOR50 use a spectrum of high frequencies, such that only a small part of the energy delivered translates in temperature rise; in fact, the main component of the energy delivered to the tissue reacts directly at the level of the molecular bonds, breaking them in a non thermal manner. This effect can be obtained thanks to the fact that the frequencies used by the RESABLATOR50 are extremely high and have values for which the resonance induced on the bonds themselves is maximum.

The incision achieved that way is very delicate because the electrode never exceeds the temperature of tissues necrosis and does not result in side laceration (the cut does not require manual pressure from the surgeon).

Instead, with the RESABLATOR50, when selecting the coagulation function, the incident quanta are characterized by an energy value slightly off-resonance. That way, since the energy associated to the Quantum is slightly different from the energy of the tissue molecular bonds, the result is a modest rise of temperature, enough to trigger the coagulation process by protein denaturation, but lower than the temperature of necrosis. This type of natural coagulation has the great advantage that it does not cause the blood vessels to collapse, as opposed to what happens in the case of traditional electric or radio-surgery scalpels.

Associating the simultaneous coagulation feature (mixed effect) achieved by causing the current to pulsate, to the cutting function, you can obtain even better results, because the small blood vessels opened by the incision are rapidly closed.

The RESABLATOR50 was built with cutting-edge electronic components and technologies, some of which have been specifically tested and in strict compliance with the international standards: the safety systems are double and are controlled by two microprocessors; the leakage currents on the patient, on the cover and to the ground are much lower than those allowed by the IEC standards.

The RESABLATOR50 uses handles in which electrodes are inserted for cutting, cutting+coagulation, haemostasis, electric shock and an insulated neutral plate to be placed in contact with the patient as close as possible to the intervention area. We can also supply disposable insulated plates, specifically designed to be used with the Vesalius® models.

Although the contact between plate and patient is not perfect, the high frequency and the features of the insulated plate allow the current to flow efficiently, avoiding burns; therefore, there is no need for the device to check the contact between plate and patient.

In the light of these considerations, the RESABLATOR50 are unique and versatile units, that allow surgical procedures on soft tissues in a precise and effective manner, while ensuring extreme delicacy and respect for the tissue itself and the neighbouring structures.

2.2.1 TECHNOLOGICAL PRINCIPLES OF THE VESALIUS® RESABLATOR50

The energy developed by the RESABLATOR50, in **RESABLATION** mode, when displayed on the oscilloscope screen, is:

- d) a modulated, high frequency sine wave (cutting + coagulation).
- e) a modulated, high frequency sine wave (coagulation).

The best results are obtained, depending on the surgeon's technique, by combining the most suitable electrode with the most suitable power output. For example, when making an incision, speed and movement are factors that allow the reduction of dehydration and tissues shrinkage.

The operating speed depends on the power set on the device, on the impedance and on the contact area of the electrode used. For example, with a large loop the current density is lower, due to the wide circumference of the electrode and results in greater energy leakage. Instead, if the loop is small, the current density is higher around the smaller interception area; therefore less power will be necessary to achieve the result desired.

When working on a "wet" tissue you have low impedance, while a dry skin tissue has high impedance.

The energy developed by the RESABLATOR50, in **RESABLATION** mode, when displayed on the oscilloscope screen, is:

- f) a symmetrical (alternated) sine wave, with zero mean value.

For further information see paragraph 4.5.1.

2.2.2 CHARACTERISTICS OF THE RESABLATOR50

The RESABLATOR50, like other VESALIUS models are made for cutting and coagulating tissues during electrosurgical procedures. Two types of accessories can be connected to the front panel: one monopolar outlet for connection to the **Resaflex™** probes (monopolar coagulation function) and one bipolar outlet (bipolar blend function) for **Resadisc™** probes.

Each function can be adjusted independently by means of potentiometers placed on the front panel.

The RESABLATOR50 is equipped with safety systems compliant with the international standards for electrosurgery devices: it has a safety circuit whose purpose is to disable the output and signal, by means of an audible and visual alarm, any connection failure with the neutral plate.

Also, it can be supplied with foot switches both electrical and pneumatic; in standard configuration it has pneumatic foot switches, which eliminate the inconvenience of having live electric cables on the floor.

The new RESABLATOR50 is also equipped with a stimulator inside, able to facilitate the identification of the position of the affected nerve before proceeding with surgery. Such purpose is achieved by means of low intensity electrical pulses, which produce a muscular reaction of the tissue adherent to the nerve, in proximity of which stimulation is applied.

2.3 TECHNIQUES AND FUNCTIONS USED

MONOPOLAR FUNCTIONS

Coagulation

The electric field generated by the QUANTUM MOLECULAR RESONANCE technology in COAG function transmits a flow of energy to the tissue, that allows the haemostasis of the blood vessels by protein denaturation of the fibrinogen, which, by turning into fibrin, closes the vessel without causing them to collapse, thereby reducing the thermal damage and the necrosis of the tissue itself.

For coagulation, usually a wide surface electrode is used and we recommend touching the tissue gently and intermittently, without putting pressure on it, until bleeding stops. It is important to keep the electrode moving, in order to avoid overheating the tissue.

BIPOLAR FUNCTIONS

Cut-coagulation

With the CUT-COAG (BLEND) function you lose a small part of the pure CUT purity to the benefit of a simultaneous coagulating effect. In this case, the effect is reduced when compared to the coagulation achieved with the "COAG" function.

When the tissue is areolar, an imperceptible coagulation film is formed, which does not interfere with the normal healing that is the primary objective and which falls spontaneously when healing is complete.

When dosing correctly the power output and using a bipolar clamp with contact surface as small as possible, the CUT-COAG function in bipolar mode allows you to achieve a delicate separation and cutting of the tissue. The CUT-COAG is also used to perform ablation or tissue lysis.

The power to be used must be adapted to the type of application so as to achieve the desired effect on the tissue, without causing sparks.

The power scale has been studied so that in its first half the variation is minimal, presenting excellent resolution; this part of the scale is perfect for very delicate tissues; in the second half it has a more significant variation, being indicated for other types of tissues.

THE OPERATING FEATURES ARE GUARANTEED ONLY WHEN USING ORIGINAL ACCESSORIES.

2.4 THE PROFILE OF THE OPERATOR AND THE ENVIRONMENT WHERE IT IS USED

The RESABLATOR50 is an active invasive device, which, by requirements and conformity, falls in the category of high frequency electrosurgery devices and as such it can only be used by qualified medical staff, trained in the particular techniques and surgical procedures to be followed.

The RESABLATOR50 must be used in a healthcare environment (operating room, medical clinic etc) compliant with the regulations on electrical systems and the environmental and electromagnetic conditions (temperature and humidity) listed below.

The positions of the operating physician and of the patient, from an operational point of view, are common to those of any electrosurgery unit. There are no entrapment areas.


The RESABLATOR50 must only be used by qualified personnel and only after they have fully read the safety and operating instructions given in this manual, supplied by the staff who marketed the device and its accessories. The medical staff, which is the end user of this system, is required to have also certain skills and knowledge of anatomy and neurophysiology and also to know the operational protocols used to perform peripheral nerve blocks.

Telea offers the possibility to attend training courses, in order to supply the user with the necessary knowledge and skills needed for the proper use of the device. Duration and frequency of the said courses are customized, based on the customer needs.


3 DESCRIPTION OF THE DEVICE

3.1 SYMBOLS

This paragraph explains the meaning of the symbols used in the device and on the labels, in accordance with the international standards.

	Fixed part, type CF, protected against the effects of defibrillator discharge (RESABLATION mode)		Keep away from sunlight
	BF-type fixed part (NERVE STIMULATION mode)		Alternating Current
	Caution		Refer to directive WEEE
	Insulated HF Patient circuit (floating)		Number in the catalogue
	Monopolar mode		Batch code
	Neutral plate		Serial number
	Bipolar mode		Rotary knob variator
	Plate detachment		Keep dry
	European Community Marking		Temperature limitation
	Equipotentiality		Use before
	Non-ionizing radiation		Do not re-use
	Fuse		Ethylene oxide sterilization
	Manufacturer		Non sterile
	Warning: read the instructions for use		Latex free
	Read the user manual		

3.2 FRONT PANELS DESCRIPTION

#	FRONT PANEL
1	IGNITION LIGHT SWITCH
7	MONOPOLAR CONFIGURATION ENABLING BUTTON
9	BIPOLAR CONFIGURATION ENABLING BUTTON
16	FAULT LED
17	PLATE DETACHMENT INDICATOR LED
19	CONNECTOR FOR MONOPOLAR HANDLE
21	CONNECTOR FOR INSULATED PLATE
22	CONNECTOR FOR BIPOLAR
23	NERVE STIMULATOR ENABLING BUTTON
24	MOTOR STIMULUS SELECTING BUTTON
25	SENSORY STIMULUS SELECTING BUTTON
26	SCREEN FOR DISPLAYING THE VOLTAGE OF THE STIMULUS IN NERVE STIMULATION MODE
27	SCREEN FOR DISPLAYING THE OUTPUT POWER (RESABLATION MODE)
28	LED INDICATOR FOR MONOPOLAR OUTLET OPERATION (RESABLATION MODE)
29	LED INDICATOR FOR BIPOLAR OUTLET OPERATION (RESABLATION MODE)
30	LED INDICATOR FOR OPERATION IN THE NERVE STIMULATION MODE
B	CONTROL KNOB FOR ADJUSTING THE VOLTAGE IN THE NERVE STIMULATION MODE
G	CONTROL KNOB FOR ADJUSTING THE VOLTAGE IN RESABLATION MODE

3.3 DESCRIPTION OF THE FRONT PANEL CONTROLS

The control panel can be divided into four sections:

- selector switches and buttons section
- outlets section;
- adjustments section;
- display section.

3.3.1 SELECTOR SWITCHES AND BUTTONS SECTION

ON/OFF switch	
	Therefore, with the ignition switch (1) in ON position, the device is supplied with power, the screens (if any), the LEDs inside the keys and the FAULT LEDs (16) light up. While these LEDs remain on, the system performs a test of the electronic circuits to ensure proper operation; in case of malfunction, only the FAULT LED (16) stays on. To restore operation you need to press the button SET (2), to perform the test again.

Function selecting buttons in NERVE STIMULATION mode

To choose one of the two stimulation frequencies you need to first enable **NERVE STIMULATION** mode, by pressing the “ENABLE NERVE STIMULATOR” (23) button. Enabling occurs when the LEDs inside the button (23) light up. In this case, **RESABLATION** mode is, consequently, disabled.

These buttons are related to the selection of the desired stimulation frequency. The buttons “2Hz MOTOR” (24) and “50Hz SENSORY” (25) cannot be selected simultaneously, therefore by enabling one frequency you disable the other:

- 2Hz MOTOR: to be used for motor nerves stimulation
- 50Hz SENSORY: to be used for sensory nerves stimulation

Enabling of a function occurs when the LED inside the button (23) light up.

These buttons are related to the selection of the desired configuration. The monopolar or bipolar configurations can be selected only in **RESABLATION** mode. Therefore, if you want to go from the monopolar NERVE STIMULATOR function to the bipolar NERVE STIMULATOR function or viceversa, you must first select and enable (by pressing the foot switch) that particular configuration in **RESABLATION** mode. The buttons MONOPOLAR (7) and BIPOLAR (9) cannot be selected simultaneously and therefore, by pressing one of the two buttons you disable the other one.

The RESABLATOR50 operates in **NERVE STIMULATION** mode or in **RESABLATION** mode alternately. In both operation modes the device uses the same outlets and the same accessories (be they active or neutral).

Function selecting buttons in RESABLATION mode

To choose one of the two power functions (**RESABLATION** mode) you need to first disable **NERVE STIMULATION** mode, by pressing the “ENABLE NERVE STIMULATOR” (23) button. Disabling occurs when the LED inside the button (23) turns off. In this case, **NERVE STIMULATION** mode is disabled and therefore you can select and use the two configurations: MONOPOLAR and BIPOLAR.

These buttons are related to the selection of the desired configuration. Enabling of a function occurs when the LED inside the button (23) lights up. The buttons MONOPOLAR (7) and BIPOLAR (9) cannot be selected simultaneously and therefore, by pressing one of the two buttons you disable the other one.




Select MONOPOLAR mode in case you intend to use a monopolar technique that requires the use of accessories connected to outlets (19) and (21). Or, select BIPOLAR mode in case you intend to use a technique that requires the use of a bipolar accessory connected to outlet (22).

The RESABLATOR50 operates in **NERVE STIMULATION** mode or in **RESABLATION** mode alternately. In both operation modes, the device uses the same outlets and the same accessories (be they active or neutral).

3.3.2 OUTLETS SECTION

The outlets available on the panel are correlated with the modes of operation available on the RESABLATOR50 device. In fact, both **NERVE STIMULATION** mode and **RESABLATION** mode can be used both with a monopolar configuration (one pole active and one neutral) and with bipolar configuration (two poles). In general, they can be divided in three sections:

- MONOPOLAR;
- NEUTRAL;
- BIPOLAR.

MONOPOLAR connector	
	The MONOPOLAR outlet can be easily identified on the front panel; it is found on the right side of the panel. Said outlet is available for the monopolar mode of operation in the function COAG (in RESABLATION mode) and for the NERVE STIMULATION mode. Therefore, you can use both modes during the same intervention without having to intervene on the accessories connected to the device. N.B.: To use the COAG function you need to connect the handle to the specific power plug (20), found on the front panel of the device. This function is only enabled for use with footswitch-controlled monopolar accessories. N.B.: Do not connect and do not position directly on the patient or next to the same other cables or accessories.
NEUTRAL connector	
	The NEUTRAL connector (21) is used to connect the disposable plates adapter and the conductive bracelet code 2505016, necessary to operate both in NERVE STIMULATION mode and in RESABLATION mode with the monopolar configuration. N.B.: Do not connect and do not place directly on the patient or next to the same other cables or accessories.
BIPOLAR connector	
	The BIPOLAR outlet is available for the bipolar mode of operation in the function CUT/COAG (in RESABLATION mode) and NERVE STIMULATION mode. Here too you can use both modes during the same intervention without having to intervene on the accessories connected to the device. N.B.: To enable the function CUT/COAG it is imperative to connect the accessory to the specific power plug (22) placed on the front panel of the device and disable mode1. Enabling is achieved by using the single foot switch supplied. N.B.: Do not connect and do not position directly on the patient or next to the same other cables or accessories.

3.3.3 FUNCTIONS SETTING SECTION

The control and setting section is found on the central part of the front panel and is divided in two sub-sections easily identifiable (with black background) referring to the two modes in which the device can operate: **NERVE STIMULATION** mode and **RESABLATION** mode.

The “ENABLE NERVE STIMULATOR” button (23) distinguishes between these two modes, that is, when this button is engaged the device operates in **NERVE STIMULATION** mode, while when it is disengaged the device is set for **RESABLATION** mode and therefore for power delivery. The status of the “ENABLE NERVE STIMULATOR” button is indicated by a LED found inside the button itself; when it is turned on (green light) **NERVE STIMULATION** mode is enabled, while when it is turned off, **RESABLATION** mode is enabled.

For the next step, namely the selection of the various stimulation or power operation functions, it is necessary to press the specific buttons, selecting between 2Hz or 50Hz stimulus in **NERVE STIMULATION** mode, or between monopolar or bipolar modes of operation in **RESABLATION** mode.

3.3.4 THE CONTROL SECTION

The control section consists of two adjustable control knobs correlated to the numbers displayed on the screens above.

NERVE STIMULATOR knob (B):

The **NERVE STIMULATOR** knob (to the left with respect to the other) is used to adjust the intensity of the stimulus in **NERVE STIMULATION** mode. Turning it clockwise increases the intensity, as shown by the magnitude scale screen-printed on the panel and by the value displayed in green colour on the corresponding numeric display. The 2-digit numeric display above the knob shows the intensity of the stimulus. This knob consists of a multi-turn rotary encoder, unlike the **RESABLATION** knob, which has an operating range within the scale screen-printed in the display.

RESABLATION knob (G):

The **RESABLATION** knob (to the right with respect to the other) must be used to adjust the power output for the monopolar COAG function and for the bipolar CUT/COAG function in **RESABLATION** mode. Turning it clockwise increases the intensity, as shown by the magnitude scale screen-printed on the panel and by the value displayed in red colour on the corresponding numeric display. The 2-digit numeric display above the knob shows the level of the output power.

3.3.5 THE DISPLAY SECTION

VOLT display (26):

Provides a numerical indication of the stimulus intensity and it varies depending on the position of the multi-turn **NERVE STIMULATION** knob from a minimum of 0.1 up to a maximum of 9.9. The value showed by this display corresponds to the value in Volt actually emitted by the device on a 500Ω charge with a maximum margin of error of ±10%.

POWER display (27):

It provides a numerical indication of the power delivered based on the position of the **RESABLATION** knob for each function selected. The value showed by this display does not correspond to the value of the actual power in Watt emitted by the device, but it is a numerical indication related to the position of the knob indicator. The display shows a minimum value of 00, up to a maximum of 70 for both functions available.

The correlation between the numerical indication and the actual power in watt can be found in the graph shown in section 4.3.1.

Blue LED monopolar COAG (28):

Indicates the power delivery in **RESABLATION** mode, in the monopolar COAG function. It flashes in case of high currents (current limits) absorbed by the power stage, situation that blocks the power delivery.

If the LED flashes simultaneously with the yellow bipolar CUT LED (29) the voltage is more than 10% of the value shown on the nameplate and you need to take prompt action on the power supply.

Yellow LED bipolar CUT (29):

Indicates the power delivery in **RESABLATION** mode, in the bipolar CUT-COAG function. It flashes in case of high currents due to a decrease in impedance. If the LED flashes simultaneously with the COAG LED (28) the voltage is more than 10% of the value shown on the nameplate and you need to take prompt action on the power supply.

FAULT LED (16):

Indicates a possible malfunction during ignition, or a system anomaly. Reports an unresponsive system in the event of overheating of the power stage. Refer to section 4.4.2 for further information.


The red LED (17):

Reports a fault in the electrical continuity of the neutral return circuit (the LED flashes simultaneously with an intermittent audible alarm). This alarm disables the outlets, therefore the device WILL NOT deliver power.

The green NERVE STIMULATION LED (30):

Reports the delivery of the stimulus energy in **NERVE STIMULATION** mode, both in the function "2Hz MOTOR" and in the function "50Hz SENSORY".

3.4 REAR PANELS DESCRIPTION

#	REAR PANEL
23	VOLUME CONTROL KNOB
24	MAIN PLATE LABEL
28	POWER PLUG WITH FUSE HOLDER
29	EQUIPOTENTIAL NODE
G	CONNECTOR (BUSH) FOR SINGLE FOOT SWITCH

Under the power socket (28) there may be a two-position selector switch that allows you to select the value of the supply power in the case of devices with dual primary transformer. Refer to section 3.4.2 for further information.

3.4.1 SERIAL NUMBER

The serial number is characterized by 9 digits which contain internal information necessary for the internal quality control procedures of Telea. The date of manufacture is nevertheless visible on another label, thus maintaining compliance with the European Directive 93/42/EEC and subsequent amendments (2007/47/EC). The date is expressed in the following format:

MM/YYYY

3.4.2 CUSTOMIZATION

The "OPTIONS" field of the nameplate label contains details of the customizations made in the specific device. Said field consists of a variable number of letters, from a minimum of zero letters to a maximum of three letters.

Codice - Code		08034055290019
REF	2501002	
OPTIONS	XXX	
Matricola - Serial number		
SN	B09010101	

The first X indicates the supply power voltage, the second X indicates the presence of electrical foot switches or air-actuated foot switches, while the third X reports the presence of the impedance control board. The meaning of the letters is explained in the following legend:

The first X		The second X		The third X	
Supply power voltage 230V	No letter	Foot switches type pneumatic	No letter	Absence of the impedance control board	No letter
Supply voltage 115V	U	Foot switches type electrical	E	Presence of the impedance control board	C
Supply voltage 100V	J				
Dual supply voltage 110-120/220-230V	B				

The standard version consists of a device powered at 230V with pneumatic foot switch without impedance control board, being therefore a version with field type "- - -".

The supply voltage is also shown on the label below the power socket to which the power cable must be connected.

N.B.: Version B (dual supply voltage 110-120/220-230V) is equipped with a dual primary transformer in order to operate over an extended range of supply voltages. On the rear panel, below the power supply socket (28) there is a two-position selector switch (115V or 230V) for the selection of the correct power supply.

We recommend to check the voltage available on the power network and then set the selector switch in the right position in order to avoid causing serious damage to the device.

The warranty does not cover breakage or malfunctions due to incorrect setting of the selector switch.

3.5 DESCRIPTION OF THE REAR PANEL CONTROLS

Connector for single switch (G):

The single pneumatic foot switch must be connected to the bushing (G) by simply pressing the tube inside the bushing. In the case of the electric foot switch, connect the metal connector to the bushing (G) and then tighten the lock ring.

Volume control knob:

The knob (23) controls the volume of the outlet activation tones. The tones are different depending on the function selected and on the mode of operation used; in fact, three different operating tones can be heard: one for **NERVE STIMULATION** mode and two for **RESABLATION** mode.

Labelling and electrical safety:

The label (24) contains the nameplate data as required by the applicable standards of the IEC 60601-1 series and compliant with the requirements on labelling set out in directive 93/42/CEE and s.m.i.. Terminal (29) is the equipotential node. Plug (28) incorporates two network protection fuses and should be connected to the power cable supplied.

3.6 STANDARD EQUIPMENT OF THE DEVICES

The following standard equipment is to be considered illustrative. Refer to the updated catalogues for a complete list of all the different parts that can be connected to RESABLATOR50:

STANDARD ACCESSORIES		
Q.ty	Description	Code
20	Disposable neutral insulated plate	2503003
1	Adapter for disposable plates with elastic band	2505016
1	Single pneumatic foot switch	2504006

Each device is supplied with the power cable, depending on the country of the buyer, and with the user manual.

4 TECHNICAL SPECIFICATIONS

4.1 INPUT TECHNICAL DATA

Supply voltage	100V or 115V or 230V or 110-120/220-230V~ 50/60Hz
Maximum absorbed power	550VA

* **Warning:** On each device can be mounted four different types of transformers that allow the connection to the power supply systems above. Refer to the voltage value indicated on the label placed below the power supply socket behind the device in order to know the appropriate value for that specific device.

The double power supply model is provided with a switch with two positions that allows the selection of the power supply value; the other models are provided with only one nominal voltage - 100V or 115V or 230V. Please use the correct power supply.

4.1.1 FUSES

The power socket (28) contains an accessible fuseholder in which are placed two fuses having the following features:

- type T 3.15A and low breaking capacity for 230V networks;
- type T 8A and low breaking capacity for 100V, 115V networks and for double-supplied devices 110-120/220-230V;

If these fuses get burned, they should be replaced with similar fuses of the same model and features (voltage, current, triggering speed and interruption capacity).

4.2 OUTPUT TECHNICAL DATA

All specifications are nominal and can be changed without prior notice. For every indicated specification should be considered a $\pm 20\%$ tolerance from the declared value, considering a temperature of the environment of 25°C and a nominal line input voltage equal to the value on the label below the power supply socket (28).

4.2.1 OUTPUT POWER

Output Frequency	Spectrum at high quantized frequencies
-------------------------	--

The following output power values correspond to the output nominal values [W] concerning the indicated nominal load [Ω] and match the indications of the identification plate placed on the rear panel of the device in **RESABLATION** mode.

NERVE STIMULATION mode does not imply power output, but a low incentive voltage having the following features.

While measuring the high frequency emissions, particular attention should be paid as to ensure precision and safety; in order to reproduce such values, the test environment should be set up properly, using the appropriate tools as specified in the relative standard.

Maximum Output power		
RESABLATION mode	COAG monopolar	70W on 400 Ω
	CUT-COAG bipolar	60W on 100 Ω
NERVE STIMULATION mode	NERVE STIMULATION	10V on 500 Ω

4.2.2 OUTPUT VOLTAGE

The maximum output voltage for each functioning mode is, on average, the following:


Maximum Output voltage		
RESABLATION mode	COAG monopolar	630 V
	CUT-COAG bipolar	900 V
NERVE STIMULATION mode	NERVE STIMULATION	10 V

Warning: The active accessories should have a nominal voltage equal to or greater than the maximum indicated output voltage. The use of original accessories ensures the functional performance of the Quantum Molecular Resonance.

4.3 OUTPUT CHARTS

4.3.1 FEATURES OF THE OUTPUT POWER

The following diagrams indicate the output power for **RESABLATION** mode depending on the adjustment taking into consideration the specific load. These are indicative curves; please consult the testing report (Mod. 10-02) for the quantitative curves.



The diagrams indicating the output power after full and half adjustment taking into consideration a load on a range for the load resistance between 10Ω and 2000Ω for monopolar outputs and between 10Ω and 1000Ω for bipolar outputs are reported within the testing report (Mod. 10-02) delivered along with the device.

4.3.2 FEATURES OF THE OUTPUT CURRENT AND VOLTAGE

The following diagrams indicate the output current and voltage in **NERVE STIMULATION** mode depending on the encoder adjustment (B) taking into consideration a characteristic load of 500Ω .

The ideal curve corresponds to the value indicated on the display (26), while the real curve corresponds to the value expressed in voltage/current actually generated by RESABLATOR50.

4.4 DEVICE FEATURES

SIZE RESABLATOR50	
Width	340mm
Height	170mm
Depth	320mm
Weight	13.5kg

Features of RESABLATOR50:

- QUANTUM MOLECULAR RESONANCE device with patented technology;
- “soft” nerve stimulation thanks to the particular wave form;
- high precision and accuracy obtained with the quartz;
- electronic components in solid state;
- adjustment of the output power using potentiometers connected to digital displays;
- adjustment of the power of the nerve stimulator using precision encoder and digital display;
- power memory set on 2 display even if off and/or not powered
- high insulation floating patient circuit with outputs protected against defibrillator use.
- double safety circuits for inner faults control, overheating and patient plate;
- self-test at every ignition;
- acoustic and visual signalling indicating power activation; different sound will be played depending on the function chosen.
- acoustic and visual signalling indicating the activation of the nerve stimulation pulses;
- system made of 3 microprocessors for controlling the system, display and a dedicated one for NERVE STIMULATION mode;
- noise network filters;
- container made of hard polystyrene, scratch-resistant and impact-resistant, with internal shielding for electromagnetic compatibility in compliance with all applicable UL 94 V-0;
- functions selectable via LED membrane buttons that indicate the selected function;
- continuous forced cooling;
- filters towards the applied part to remove any low frequency components, that may cause neuromuscular stimulation;
- pneumatic pedal; electrical pedals protected against entry of dust and liquids; optional;
- it complies with the applicable standards in force regarding fundamental safety and essential performance of the electro-medical devices series IEC EN 60601-1 (IEC EN 60601-1-2, IEC EN 60601-2-2, etc);
- during the tests implied by the previously mentioned standards, the device has been subjected to the most unfavorable conditions, both from the climatic point of view (see maximum allowed temperatures and humidity levels) and from the supply point of view (for example 100V supply voltage for absorption measuring and 60Hz frequency to measure the dispersion currents, etc.).
- in compliance with European Directive 93/42/CEE and smi (2007/47/CE). Medical device class IIb.

4.4.1 SELF-TEST AT START-UP

At each start-up, the RESABLATOR50 tests its electric circuits as to ensure proper operation and simulates faults as to check the correct intervention of the control systems. During this testing period, the FAULT LEDS and the digital displays are on; this way the operator may check the test time and the proper functioning of the LED. The LED CUT, COAG and plate detach may also be controlled by powering them up, first with the plate connected and then powering up the RESABLATOR50 without the plate being connected.

4.4.2 ALARMS

Temperature check:

The device is provided with inner probes for temperature check.

An intensive use of the device, without respecting the declared ON/OFF cycles and the data on the identification plate, may lead to a dangerous increase in the internal temperature. The intervention of the control devices deactivates the outputs and on the display will appear the fault condition due to the temperature (FAULT LED and all segments of display LED on).

Therefore, you should wait for the RESABLATOR50 to cool down, and after about 10 minutes you may turn it off and then restart it. In fact, the device should be let on to speed up the process of lowering the internal temperature.

Neutral plate check:

If the neutral plate is not connected, when the monopolar inputs are activated, the plate detach LED will turn on (17) accompanied by an intermittent sound. The alarm continues until the plate is connected.

Current Control:

Work situations with low-or very low impedance or in case of parasitic capacitances can lead to an additional increase of the current in the power circuit at high frequency.

The intermittent operation of the blue monopolar LED COAG (28) or of the yellow bipolar LED CUT-COAG (29) indicates the intervention of the current control and deactivates the outputs in intermittent mode. This signal will automatically deactivate as soon as the normal operating conditions are restored. If the anomaly persists, the power should be decreased until the LED generates a steady light.

Voltage check:

The voltage check is performed as to make sure that the value of the output powers is kept within the allowable limits. If the voltage increases beyond the value of +10%, a visual and acoustic alarm will be triggered (simultaneous flashing of LED (28) and (29)) and the device will be shut-down. It is recommended the use of a stable voltage value.

4.4.3 PROTECTION FROM DEFIBRILLATION

RESABLATOR50 is provided with applied parts, protected against the discharges that take place during defibrillation. After 10 seconds from the exposure to the defibrillation voltage, the device returns to comply with the requirements of the standards regarding safety and essential performance.

4.5 NERVE STIMULATOR TECHNICAL DESCRIPTION

4.5.1 GENERAL INFORMATION

Below are indicated general informations regarding the output signal generated by RESABLATOR50 when set to **NERVE STIMULATION** mode:

Frequency of the impulse: There are two stimulation frequencies available: 2Hz or 50Hz that can be selected using the two buttons (2Hz MOTOR" (24) and "50Hz SENSORY" (25)) on the front panel (see par. 3.3.1).

Duration of the impulse: The signal impulse is made of an alternate sinusoidal wave; the duration of the semi-wave is 1 ms, therefore the overall duration of the impulse is 2 ms.

Amplitude of the impulse: The signal amplitude can be adjusted up to a maximum intensity of 20mA (that correspond therefore to 10V on 500Ω). The maximum value of 20mA is never exceeded, even if the applied load varies.

Display: The amplitude of the impulse can be viewed on two digital displays with 7 segments. The minimum value displayed is 0.1 up to a maximum of 9.9 with a pitch (resolution) of 0.1 in relation with the position of the encoder (B). This indication concerns a voltage that corresponds to a characteristic charge of 500Ω. Given the variable impedance of biological tissues, the displayed value should be considered as an indication of intensity rather than an exact value in volts of the output signal according to the graphs shown in the next paragraph.

Adjustment: The amplitude of the signal can be adjusted from a minimum value of 0.1 up to a maximum value of 9.9 with a resolution of 0.1 corresponding to 100 adjustment pitches. The display shows the voltage corresponding to a charge of 500Ω. For charges other than 500Ω the adaptive system of the NERVE STIMULATOR still allows the adjustment of the impulse amplitude safely, because the maximum intensity set (value of 9.9 appears) does not exceed 20mA and decreases according to the encoder adjustment (B).

N.B.: Nerve stimulation may occur in monopolar or bipolar mode, depending on the last function used (with power supply) in **RESABLATION** mode.

4.5.2 WAVEFORM


The waveform of the stimulation signal is a sinusoidal impulse with a duration of 1 ms for the positive half-wave and 1 ms for the negative half-wave, for a total of 2 ms. After the impulse comes an off phase, that lasts for a period of time that varies depending on other settings made by the user.

The off phase has a duration of:

- 20-2=18 ms if the 50Hz function is set (duration of the waveform 20 ms)
- 500-2=498 ms if the 2Hz function is set (duration of the waveform 500 ms)

An example of representation of the stimulation signal is shown in 4.5.2-1. It refers to a charge of 500Ω.

The "smooth" sinusoidal waveform allows a more gentle stimulation of the tissue, compared to a squared waveform, allowing a greater stimulation intensity, facilitating the diagnostic.


4.5.2-1

4.5.3 SAFETY AND PRECISION

Thanks to its adaptive system, **NERVE STIMULATION** mode allows the stimulus adjustment even for charges different from 500Ω.

As described in the previous paragraph, RESABLATOR50 performs an automatic control of the output impedance and limits the current (stimulus intensity) to the maximum value of 20mA into the charge area from 0 to 500Ω and then decreases gradually while the impedance increases beyond 500Ω.

Figure 4.5.3-1 shows the average performance of the stimulus current on the characteristic charge of 500Ω with the encoder (B) adjusted at maximum and at half of the maximum value.

4.5.3-1

4.6 CLASSIFICATION

The device is classified as:

- Class I device with type “CF” applied part when used in **RESABLATION** mode;
- Class I device with type “BF” applied part when used in **NERVE STIMULATION** mode;
- device with casing, not protected against the ingress of water (IPX0);
- the device is not suitable for use in the presence of a flammable anaesthetic mixture with air, oxygen or nitrous oxide;
- device for intermittent (non-continuous) operation, with ON and OFF cycles; be sure to read the ranges reported on the label placed on the back of the unit.

4.7 INTENDED USE

The device is intended, but not limited to the following applications:

- Endoscopic epiduroscopy interventions
- Ablation and lysis of adhesion tissue

4.8 ELECTROMAGNETIC COMPATIBILITY

All design solutions for RESABLATOR50 led to the development of medical devices able to be able to function properly in a controlled and appropriate electromagnetic environment, in compliance with minimum safety distances from mobile terminals (as indicated by standard EN 60601-1-2).

The user must therefore make sure that the electromagnetic environment has the features described below:

- c) Electromagnetic emissions:

Emission Test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Unit 1	RESABLATOR50 uses RF power only for internal functions. As a consequence, the RF emissions are very low and do not disturb the nearby electronic devices.
RF emissions CISPR 11	Class B	RESABLATOR50 can be used in any environment, including domestic environments and environments powered directly from a low voltage public network for domestic use.
Harmonic emissions IEC 61000-3-2	Class A	
Flicker voltage/emissions variations IEC 61000-3-3	compliant	

The presence of other cables or active accessories in the same environment, may generate disturbances that can be observed during the operation of RESABLATOR50.

- d) Electromagnetic immunity:

Immunity test	IEC 60601-1-2 Level test	Level conformity	Electromagnetic environment
Electrostatic discharges IEC 61000-4-2	6kV contact discharge 8kV air discharge	IEC 60601-1-2 Level test	Floors should be made of wood, concrete or ceramic. If the floors are made of synthetic material the relative humidity should be at least 30%
Burst/quick electric transients IEC 61000-4-4	2kV for power supply networks 1kV for input/output lines > 3m	IEC 60601-1-2 Level test	Mains power quality should be that of a typical commercial or hospital environment.
Over-voltage IEC 61000-4-5	1kV of differential mode 2kV of normal mode	IEC 60601-1-2 Level test	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_n for 0.5 cycles 40% U_n for 5 cycles 70% U_n for 25 cycles 0% U_n for 5 s	IEC 60601-1-2 Level test	Mains power quality should be that of a typical commercial or hospital environment. If the user wants Vesalius® to remain active even in case of power outages in the network, then connect it to an uninterruptible power supply or a battery.

Immunity test	IEC 60601-1-2 Level test	Level conformity	Electromagnetic environment
Magnetic field at the mains frequency (50/60Hz) IEC 61000-4-8	3 A/m	IEC 60601-1-2 Level test	The levels of the magnetic fields of the network frequency must be those typical for a commercial or hospital environment.
			The portable and movable RF communication devices must be used at a minimum distance from all parts of Vesalius®, including cables. Such distance is calculated based on the equation concerning the transmitter frequency.
RF irradiated EN 61000-4-3	3 V/m 80 MHz to 2.5GHz	3 V/m	Recommended separation distance: $d = 1.2 \times \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5GHz
RF conducted EN 61000-4-6	3 V 150 KHz to 80MHz	3 V rms	Recommended separation distance: $d = 1.2 \times \sqrt{P}$
			Where P is the maximum output power of transmitter (W) based on the specifications provided by the manufacturer of the transmitter and d is the recommended separation distance (m). The field forces for fixed RF transmitter, determined as a result of a local electromagnetic research, must be lower than the conformity level for each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

RESABLATOR50 should be used in an electromagnetic environment in which the irradiated RF disturbances may be controlled.

The distance between the mobile terminals or the wireless technologies and the medical device is one of the factors that increase the malfunctions risk. The user may avoid the electromagnetic interferences by maintaining a minimum distance between the RF portable and mobile transmitters and RESABLATOR50. This distance is calculated based on the maximum output power of the communication device.

Maximum nominal output power of the transmitter [W]	Separation Distance According to Transmitter Frequency [m]		
	150 kHz to 80 MHz $d = 1.2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \times \sqrt{P}$	800 MHz to 2.5GHz $d = 2.3 \times \sqrt{P}$
0.01;	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23;

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated from the equation applicable to the frequency of the transmitter of interest, where P is the maximum nominal output power of the transmitter in watts (W), as indicated by the manufacturer of the transmitter.

Notes:

- (1) at 80 MHz and 800 MHz, will be applied the separation distance corresponding to the higher frequency range.
- (2) These guide lines may not be applicable in all situations. Electromagnetic propagation is affected by absorption and it is a reflection phenomenon corresponding to structures, objects and people.

4.9 FEATURES OF ACCESSORIES

For a full list of accessories, additional equipment (e.g. Smoke-Buster fume extractor, trolleys), please refer to the catalogue. Each code corresponds to a different accessory; refer to the item number shown on the packaging of individual accessories.

Warning: Maximum nominal voltage for all accessories is 750V; it is recommended to use the accessories adjusted at maximum 80% of the available power. Using the accessories at maximum output power may lead to a high frequency output voltage greater than the nominal voltage of the accessories that may damage the same.

Due to the convenience of technology and frequencies used, the materials of which the active separable electrodes are made (both disposable and reusable) play an important role in exploiting the QUANTUM MOLECULAR RESONANCE phenomenon and in achieving the desired performance, even concerning the safety of the patient.

The mechanical fall test provided by the general standard for electrical safety and essential performance have been carried out from a height of 1.5 meters worsening the condition of normal use specified in the regulatory requirement.

Telea ensures compliance and correspondence with the technical specifications described in this manual only if using original accessories, designed for exclusive use with the RESABLATOR50.

5 INSTALLATION, USE AND WARNINGS

5.1 CONDITIONS OF USE, TRANSPORT AND STORAGE

The conditions of use and storage of the accessories may vary from those indicated below, that refer to the device. Please consult the label and the instructions inside the packages of the original accessories.

	ROOM TEMPERATURE	RELATIVE HUMIDITY	ATMOSPHERIC PRESSURE
CONDITIONS OF USE	+10°C (50°F) ÷ +40°C (104°F)	30% ÷ 75% not condensed	700hPa ÷ 1060hPa
TRANSPORT AND STORAGE	0°C (32°F) ÷ +50°C (122°F)	20% ÷ 80% not condensed	500hPa ÷ 1060hPa

5.2 INSTALLATION OF THE DEVICE

During RESABLATOR50 installation, please comply with the provisions below.

- Place the device on a solid surface, not bulky, away from the patient and from areas potentially subject to accidental bumps so as to avoid crushing or falling risks.
- For proper ventilation of the device, it must be far enough away from the walls that obstruct air circulation and free from any encumbrances or objects placed on the bottom or side that can facilitate the heating of the internal circuitry.
- Arrange the equipment so that the power cord can be easily disconnected from the power plug in case of need for isolation from the mains.

5.3 CONNECTION OF ACCESSORIES

Single pedal:

The pedal connection cable must be inserted in connector (G) (black coloured) placed on the rear panel of the device.

Once the pedal is connected to the machine, select the desired mode using the appropriate key (23), then choose the type of function you want and activate it by exerting pressure on the pedal.

Active accessories:

The active accessories should be connected to the front panel of the device. After connecting the panel, insert in the monopolar output (19) of the front panel the monopolar accessory to be used (e.g. RESAFLEX probe). Similarly, if using bipolar techniques, insert the bipolar accessory to be used in the bipolar output (21) of the front panel.

Insulated neutral disposable plate:

Insert the adapter of the plate into the NEUTRAL output of the front panel (21) and lower the lever in order to attach the metal plate permanently. Please see the next paragraph for more information.

Elastic bracelet:

A conductive elastic bracelet is supplied with the connector for the insulated neutral plate and allows nerve stimulation. The elastic strap is detachable from the adapter because connected via a banana plug.

5.4 GENERAL INSTRUCTIONS OF USE FOR THE NEUTRAL PLATES

RESABLATOR50 use disposable neutral plates with high level of insulation, that reduce the risk of burns.

5.4.1 DISPOSABLE PLATE

Apply the plate in a very vascularized point and as close as possible to the area of surgery. The current pathway within the body of the patient should be as short as possible and should proceed diagonally, without flowing in the transverse direction through the chest. Avoid, if possible, areas with particularly thick hair otherwise shave the area. Bone prominences, scars, metal implants, electrodes of other equipment (for example, monitoring electrodes, ECG, EEG, etc.) must be absolutely avoided.

Clean and dry the area of application; when using disinfectants make sure not to dampen the adhesive glue. Residual moisture must be dried quickly and completely.

Apply the plate starting from the side of the width and continuing toward the tongue (length), pressing with your fingers the sticky edges, massaging then the entire area in contact with the patient. Make sure that the entire surface has a stable contact with the skin.

Place the adapter for disposable plates so that it does not touch the patient. The cable must not be wrapped and must not come into contact with the patient or other cables.

At the end of the intervention, gently remove the plate in the direction opposite to that of the application. Once applied, do not replace the plate in the other position but use a new neutral plate.

WARNING: If during the intervention the operator has the feeling that the cutting and / or coagulation are no longer the initial ones, IMMEDIATELY stop and do not resume the operation without carefully checking that the conductive surface is entirely in perfect contact with the patient's skin. Incorrect positioning of the neutral plate may cause burns or poor electrical performance of the generator.

5.4.2 ELASTIC BAND

The conductive elastic bracelet is not sterile when supplied. It allows stimulation of the nerve of interest using the same bipolar accessory used in **RESABLATION** mode. It is of immediate and simple application and has a connector button to allow its positioning without the cable. After positioning the bracelet, connect it to the machine using the provided cable.

The bracelet can be opened and can be placed as to be completely in contact with the patient at the point of application. It is recommended not to apply the elastic band on the skin if it presents abrasions, rashes or sores and also to avoid the band critical points for the placement of the disposable plate.

The bracelet should be placed on the limb closest to the area to be treated.

5.5 PREPARATION FOR USE

5.5.1 USE AND START-UP PROCEDURE

Given the double functionality of the RESABLATOR50, the procedure specified below applies for both operational modes.

Before starting the device, check that the voltage indicated on the label on the rear panel matches the mains voltage and that the switch position is correct (for devices with dual power supply). A stabilised voltage should be used.

- 11) connect the power supply cable to the socket (28) on the rear panel of RESABLATOR50 and to the mains;
- 12) connect the pedal actuator bushing (G) on the rear panel, marked with "FOOTSWITCH socket";
- 13) connect the monopolar accessory to the MONOPOLAR connector; the insertion of the accessory is facilitated by the centring sign on the front panel;
- 14) connect to the connector (21) the bracelet disposable adapter plates for the use of the monopolar function and of **NERVE STIMULATION** mode in monopolar configuration. For bipolar function, connect the bipolar accessory to the connector (22);
- 15) connect the conductive elastic bracelet to the patient, as close as possible to the treated area, as reported in the specific paragraph. If the bracelet is not connected, please connect the parts;
- 16) turn on RESABLATOR50 bringing the switch (1) on "ON";
- 17) choose the desired operating mode by activating or deactivating the button (23), choosing therefore between **NERVE STIMULATION** mode (button (23) activated) or **RESABLATION** mode (button (23) deactivated).
- 18) then choose the desired stimulation frequency, 2Hz (24) or 50Hz (25) (blue buttons) for **NERVE STIMULATION** mode or select buttons (7) or (9) for **RESABLATION** mode selecting between monopolar or bipolar function.

The monopolar or bipolar configuration can be selected only in **RESABLATION mode. Therefore, if you want it to switch from monopolar NERVE STIMULATOR to bipolar NERVE STIMULATOR or vice versa, that configuration must first be selected and activated (pedal pressure) in **RESABLATION** mode.**

- 19) adjust the power level / intensity of a stimulus via the two VOLT levers (lever (B) for **NERVE STIMULATION** mode) and POWER lever (G) for the **RESABLATION** mode). We recommend starting with a low value and increase it gradually following the relative indications.
- 20) Activate power (**RESABLATION** mode) or stimulation (**NERVE STIMULATION** mode) using the pedal actuator;

N.B: Nerve stimulation may occur in monopolar or bipolar channel, depending on the last function used (with power supply) in **RESABLATION mode.**

If using **RESABLATION mode with bipolar accessory, one pole of the accessory will be the active pole and the energy of the stimulus will always close on the elastic bracelet that acts as a return electrode.**

5.5.2 SHUTDOWN PROCEDURE

RESABLATOR50 can be shut down by bringing switch (1) on "OFF".

When used for subsequent interventions within a short period of time, it is recommended not to turn off the device but leave it on stand-by to allow complete cooling of the same.

After use store with extreme care all engaged reusable accessories and proceed to the disposal of the disposable ones, in compliance with the standards in force.

5.6 SAFETY PRECAUTIONS

In addition to the security measures that require the intervention of the alarm systems described in paragraph 4.4.2 the RESABLATOR50 provides additional security measures described below.

Control of Output Parameters

In order to avoid accidental activation or selection, it is not possible to change the settings if the pedal is pressed and the outputs are energized. In **NERVE STIMULATION** mode the value of the stimulus can not be changed during the delivery (pedal pressure) while it is possible to vary the output power while generating power in **RESABLATION** mode.

Obviously, the supply is interrupted as soon as you take your foot off the pedal actuator.

6 CLEANING, DISINFECTION AND STERILISATION

6.1 DEVICE CLEANING

An appropriate cleaning and maintenance, respecting the prescriptions indicated at the following points will ensure correct functioning and lasting life of the surgical instrumentation.

- Disconnect the plug from the mains before cleaning!
- Do not use flammable or explosive substances, plastic solvents or abrasives.
- The device can be cleaned using a wet cloth (for example benzalkonium chloride) and then with a dry one.
- In the case of sterilization of the environment cover the unit with a plastic sheet.
- Do not spill liquids directly on the device, do not autoclave or gas sterilize.

The device should be cleaned using a soft cloth moistened with purified or sterile water. Only detergents with neutral pH may be used. Do not spray substances against the casing and do not wet it in any way. **During cleaning, do not use excessive amounts of fluids and make sure that they do not enter into the device.** Completely dry any water residue

6.2 CLEANING AND STERILISATION OF THE ACCESSORIES

Do not use damaged tools and do not attempt to repair them by yourself. Check the integrity of the accessories before using them.

6.2.1 CLEANING

Pedal

All pedals can be washed and disinfected using gauzes soaked with disinfectant (e.g. benzalkonium chloride). There are no electrical problems in the case of infiltration of liquids; the pedals are protected against the ingress of dust and liquids.

Do not use acids or corrosive substances.

The functioning of pneumatic pedals can be checked by placing the mouth on the pneumatic valve output to check for proper air release.

Elastic band

The conductive strip can be reused and therefore it must be washed with cold water and neutral detergent after every use. Detergent solutions (non ionic surfactants) and disinfectant solutions (components based on quaternary ammonium 0.1-0.2%). Follow the instructions when using such products. After cleaning, rinse with clean water and allow drying completely before the next use.

The lifetime of this product depends on the number of uses and state of conservation and maintenance of the product. As soon as you observe deterioration of the product must be discarded.

Adapter for disposable plates

Adapter for disposable plates (Ref. 2505016) resistant to autoclaving but usually it is simply cleaned and disinfected with a cloth dipped in a special solution. Follow the instructions when using such products. Pay attention to completely dry after manual washing (to avoid the persistence of fluid in the connectors). Do not twist the cord too tightly to avoid any damage. Check the cables insulation before and after very use.

The other accessories RESABLATOR50 should be immediately disposed of after use in accordance with the procedures for the disposal of hospital waste in compliance with current regulations.

6.2.2 STERILISATION

There are no sterilization operations and reconditioning for active accessories of RESABLATOR50 because they are all disposable and non-reusable by definition.

The accessories marked as disposable are sold sterile and must not be reused as they may cause infections.

6.2.3 STORAGE

Store the accessories in a dry place, clean and free of dusts, at medium temperatures within 5°C and 40°C at a relative humidity within 20% and 80%.

7 MAINTENANCE

RESABLATOR50 does not need any particular maintenance operations; the standard regulations for electronic devices should be followed:

- periodically check the supply cable of the device, the casing and the proper operation of the front panel. Periodically check all accessories, such as active electrodes and connection cables. If they show signs of deterioration, wear, or have exceeded the number of sterilization cycles, they should be replaced;
- it is recommended to periodically check the state of the reusable insulated plate and to replace it if it shows signs of wear, uncovered metal parts or damaged integrity of white epoxy coating;
- Telea provides three forms of programmed assistance according to the needs of customers, covering the ten-year expected lifetime. After the tenth year from the date of purchase the manufacturer will offer customized solutions for assistance. The instrument can be opened only by authorised technical staff after disconnecting it from the mains.
- authorised staff should perform at least once a year a periodic inspection of the equipment; Telea recommends annual control of leakage currents at low frequency;
- Telea provides to authorised assistance centres the instructions for inspection, repair, replacement and calibration of the electronic boards and anything else requiring a specialized technical assistance;
- at the end of the operating cycle, for the disposal of the device, please refer to the hospital waste disposal procedures (refer to the indications on the last page of this manual);
- for second level maintenance of the devices of the RESABLATOR50, please consult the **VESALIUS® Service Manual**.
- please keep the original packaging materials for eventual shipping to the manufacturer or to the technical office.

For all the operations of maintenance, calibration and controls please refer to Telea or its authorized points.

8 TROUBLESHOOTING

ANOMALY	POSSIBLE CAUSE	WHAT TO DO
The unit does not start	The supply cable is incorrectly inserted or the supply voltage is missing	Check the installation of the power cable from both sides and check the presence of voltage.
	Fuses faults	To check and if necessary replace the fuse, disconnect the power cord from the mains, then remove the tray closed element combined network entry and remove the fuses. Replace the fuses with similar ones, of the same type and value
	Start button	Check that the start button R is set on ON
After start-up and self testing, the LED FAULT (16) turns on	Malfunction of the control systems	Press the button ET(2) to reactivate the test. If the problem recurs, turn the instrument off and on again after at least 10 seconds. If the problem persists, contact the Customer Service (service@vesalius.it).
LED FAULT (16) turns on after a long period of device activation.	Intervention of the temperature control systems	Wait for the tool to cool down and then press button SET (2). If that's not enough, wait to make sure that the instrument is sufficiently cooled. The cooling is facilitated if the device is left on. If the problem persists, contact the Technical Service.
By activating the device in monopolar mode LED "detachment plate" (17) turns on and a sound is heard intermittently.	Plate cable fault	Check that the neutral plate is properly connected to the tool. If the problem persists, the cable of the plate may be internally damaged; the entire piece should be replaced

9 WARRANTY

The warranty is valid for 24 months from the date of sale of the unit and covers all manufacturing defects; the warranty does not cover the accessories.

The warranty period starts from the date of purchase of the device.

The warranty does not cover:

- damages caused by improper use, negligence, normal wear and tear;
- failures due to natural disasters or fires;
- damage due to improper repairs or other work that is not authorized by the manufacturer

The warranty terms specified can not be applied to products, or part of them, repaired or modified by third parties outside the Telea factory, nor to products subjected to misuse, negligence or accidents. The shipping costs for the repair of the defective device should be borne by the buyer.

Warning: The manufacturer denies all liability for any direct or indirect damages to persons or property arising from improper use of the device or when the device has been tampered with.

Warning: The manufacturer is not liable for any damage or malfunction caused by the use of electrodes or non-original accessories.

Telea is considered responsible with regard to basic safety, reliability and performance of the electro-medical device only if:

- RESABLATOR50 was not modified, tampered with and is still sealed;
- assembly operations, fine adjustments, modifications or repairs are carried out by properly trained and authorized staff;
- the installation and preparation procedures described in this manual are followed;
- the electrical installation of the premises complies with local regulations and regulatory requirements, and
- the RESABLATOR50 device is used in compliance with the Telea instructions for use.

Please note that breaking the seal on the back or front of the equipment by UNAUTHORISED staff voids the warranty and cancels the LIABILITIES implied by the CE certification.

IMPORTANT: Please use only ORIGINAL electrodes and accessories, specifically designed for conductive characteristics and mechanical properties.

N.B.: For continuous improvements, Telea reserves the right to make changes without notice, considering valid only the latest revision of each document.

VESALIUS® is a registered trademark of **Telea**.

Telea Electronic Engineering s.r.l.

Via L. Da Vinci 13,

36066 SANDRIGO (VI)

(ITALY)

Waste Electrical and Electronic Equipment

Directives 2002/95/EC on the restriction of the use of certain hazardous substances in electrical and electronic equipment and 2002/96/EC on waste electrical and electronic equipment are designed to tackle the fast increasing waste stream of electrical and electronic equipment and complements European Union measures on landfill and incineration of waste. Increased recycling of electrical and electronic equipment will limit the total quantity of waste going to final disposal. Producers will be responsible for taking back and recycling electrical and electronic equipment. This will provide incentives to design electrical and electronic equipment in an environmentally more efficient way, which takes waste management aspects fully into account. Consumers will be able to return their equipment free of charge. In order to prevent the generation of hazardous waste, Directive 2002/95/EC requires the substitution of various heavy metals (lead, mercury, cadmium, and hexavalent chromium) and brominated flame retardants (polybrominated biphenyls (PBB) or polybrominated biphenyl ethers (PBDE)) in new electrical and electronic equipment put on the market from 1 July 2006.

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