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The **Resaflex** is a flexible monopolar electrode it consists of a wire having a length of 560 mm covered, with an insulating sheath and terminates at one end with a sphere of 0,80 mm (2,4 Fr). It is specifically designed for Epiduroscopy.

The electrode, must be inserted into the video-guided catheter **Resascope** (Please refer to the Instruction for Use relating to that device).

The **Resaflex** must be connected ONLY to the **Resablator 50** Molecular Resonance Generator.

The **Resablator 50** Molecular Resonance Generator generates a high frequency monopolar current, with a spectrum that ranges between 4 to 16 MHz that is transferred to the tissue by the **Resaflex** tip (sphere with diameter 0,80 mm (2,4 Fr) made of stainless steel).

The **Resaflex** is individually packed, sterile, single use.

Verify the integrity of the packaging before use.

If the packaging is broken or damaged the sterility of the device is not guaranteed.

PROCEDURE OF USE

- Remove the electrode from the sterile package and connect the return neutral plate to the patient.
- After entering the optic fiber into the "Optic" channel of the videoguide **Resascope** and following the instructions for use of the same, introduce the **Resascope** into the epidural space and start the investigation operating the levers *UP/DOWN* and *LEFT/RIGHT*.
- After identifying pathologies scarring or adhesions, connect the electrode to the generator **Resablator 50** using the appropriate connector (see Fig. 1). Introduce the electrode into the channel "Tool1" or "Tool2" of the videoguide **Resascope**.

Warning: The **Resaflex** is provided with a notch mark along the rigid rod adjacent to the handle, so as to signal the escape of the sphere of the **Resaflex** from the videoguide **Resascope**, when this is at the purple valve which delimits proximally the working channel.

Warning: The notch mark has purely indicative purposes and always needs to find correspondence with the endoscopic images visible on the screen. The exposure of the spherical tip of the **Resaflex** can vary by a few millimeters depending on the tightening of the purple valve located at the entrance of the **Resascope** channel.

- After checking the actual spill of the Resaflex, operate according to the depth and rotate it on itself. The objective is to achieve a direct contact between the active part (Resaflex ball) and the tissue to be removed. To perform these operations hold the probe from the black handle and help the movements also using the levers of also the videoguide.
- Only when you are assured of an effective contact between the distal sphere and the tissue to be treated, set the power of the generator using the knob of the session "Resablation-Power".
- Set an initial power between 25 W and 35 W and, after have pressed the "Monopolar" button (blue led inside the button will lightening on), start the treatment by operating the foot control pedal.
- Gradually increase the power until the desired effect and carry out the intervention (see instruction for use of the **Resascope**).

Warning: The value of working power indicatively suitable, is the one just below the value which causes the formation of bubbles (cavitation) within the epidural space.

Warning: check compatibility of the endoscopic equipment before each use.

Warning: Use only endoscopic devices classified as type CF applied part.

Warning: when energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive.

AMS Group S.r.l. s.u. is not responsible for any damage to the patient and/or medical staff caused by the anomalous or improper use of the system, wrong maneuvering, and use by unqualified staff.

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Telea Electronic Engineering S.r.l.