



**THE INNOVATIVE SOLUTION
FOR THE DECOMPRESSION
OF CERVICAL DISCS**

The Resadisc® device allows the treatment of degenerative pathologies of the cervical spine, such as:

- Disc protusions
- Small hernias
- Radicular pain caused by nerve compression/irritation

The Resadisc® allows the treatment of degenerative pathologies of the cervical spine with:


- Patented technology: **Quantic Molecular Resonance® - QMR®¹**
- **QMR®¹** with high frequency and controlled temperature (<50°C)¹
- Minimally invasive approach: **19G** introducer needle
- Bipolar device: decompression effect focussed 1mm from the electrode
- Active tip with 15° angle



CERVICAL RESADISC® Device code: 2615018

The Resadisc® kit for cervical discal decompression is composed of:

- Bipolar Resadisc® electrode for cervical discal decompression
- 19G introducer needle



RESABLATOR SMART Generator code: 2501033

Dedicated generator with Resadisc® automatic recognition and LCD display

